

medi-inn<sup>®</sup>  
CLASSIC

# medi-inn CLASSIC Protective Gowns

Professionally Dressed.  
Reliably Protected.

Optimized  
Gown Range



## Reliable Protection in Hygiene Critical Environments

Medical protective gowns must meet the highest standards: They need to reliably protect against contamination, provide comfort during long hours of wear, and allow maximum freedom of movement. Our carefully curated assortment includes visitor gowns, protective gowns with partial or full lamination, as well as isolation gowns PPE Category III – offering the right solutions for wide range of applications in hospitals, nursing care, industry, and laboratories.

## medi-inn CLASSIC – Quality You Can Trust

The medi-inn CLASSIC protective gowns stand for tested quality, functional design details, and excellent price-performance ratio. With a wide range of styles and sizes, our gowns provide optimum freedom of movement. Thoughtfully selected materials, practical closure solutions, and ergonomic cuts provide reliable support for healthcare professionals in their daily work.

## About BODY PRODUCTS

BODY PRODUCTS offers a carefully selected portfolio of more than 2,000 medical devices and healthcare items. Our range covers protection, hygiene, diagnostics, and continence care – including high-quality aids, supplies for wards and practices, as well as first aid. Another focus lies on disposable protective clothing. Innovative solutions for highly specialized medical fields such as endoscopy, surgery, and urology complete our portfolio.

## Proven Quality


As a supplier of medical devices, protective medical clothing and hygiene products, we take our responsibility seriously. That’s why we select our portfolio with the utmost care and ensure that all our products are tested for quality and performance. Our quality management system is certified according to DIN EN ISO 13485, meeting the strict requirements for medical devices. Sustainability is also a core part of our corporate philosophy: we are certified according to the environmental management standard DIN EN ISO 14001, our sustainability performance is regularly assessed by EcoVadis, and we have documented our commitment in the German Sustainability Code (DNK).



The  
SUSTAINABILITY  
Code  
Signatory 2022



# medi-inn CLASSIC Protective Gowns – Special Features



Back Binding at the Neck



Ties at the Back



Elastic Cuffs



Ties at the Side



Knitted Cuffs



Hospital



Medical Practice



Dental Practice




Laboratory



Geriatric & Nursing Care

## medi-inn CLASSIC PP Visitor Gown

The medi-inn CLASSIC PP Visitor Gown made of skin-friendly polypropylene non-woven fabric with knitted cuffs is breathable and is used as a disposable visitor gown for hygiene prevention for visitors and nursing staff or by the food industry.



Ties at the back of the neck and waist



Knitted cuffs on sleeves



- PP spunbond 16 g/m<sup>2</sup>
- Breathable
- Ties at the back of the neck and waist
- Knitted cuffs on sleeves
- Latex-free
- Category I personal protective equipment \*
- Class I medical device \*\*
- 10 pieces per bag

NEW  
Coming  
Soon



Item Number	Product	Available Size	Colour	Material & Grammage	Sales Unit (SU)
N640676	medi-inn CLASSIC PP Visitor Gown	115 x 137 cm	Light Green	PP Spunbond 16 g/m <sup>2</sup>	10 x 10 Pieces

\* According to Regulation (EU) 2016/425 \*\* According to Regulation (EU) 2017/745

## medi-inn CLASSIC Protective Gown with PE Part-Coating

The disposable gown with knitted cuffs can be used as a protective gown in basic care in the context of MRSA, MRGN and ESBL, noro and influenza infections, as well as for laboratory work and in the food industry.



**Tested to:**  
EN 14126:2003/AC:2004, ISO 16603, ISO 16604,  
EN ISO 22610, ISO 22611, EN ISO 22612

- PP spunbond with liquid-proof PE part-coating on sleeves and front
- Waist ties fasten at the sides
- Ties at the back of the neck and waist
- Knitted cuffs on sleeves
- Category I personal protective equipment \*
- Class I medical device \*\*
- Breathable and latex-free
- 10 pieces per bag

New  
Dimensions  
120 x 136 cm  
Coming Soon



Item Number	Product	Available Size	Colour	Material & Grammage	Sales Unit (SU)
N18344-y	medi-inn CLASSIC Protective Gown (with PE part-coating)	120 x 136 cm	Yellow	PP Spunbond with PE part-coating 26 g/m <sup>2</sup>	10 x 10 Pieces
N11333-y	medi-inn CLASSIC Protective Gown (with PE part-coating)	139 x 139 cm	Yellow	PP Spunbond with PE part-coating 31 g/m <sup>2</sup>	10 x 10 Pieces
N11333-y-XL	medi-inn CLASSIC Protective Gown (with PE part-coating)	140 x 160 cm	Yellow	PP Spunbond with PE part-coating 31 g/m <sup>2</sup>	10 x 10 Pieces

\* According to Regulation (EU) 2016/425 \*\* According to Regulation (EU) 2017/745

## medi-inn CLASSIC Isolation Gown with PE Part-Coating PPE Category III

The medi-inn CLASSIC Isolation Gown SELECT, made of skin-friendly polypropylene non-woven fabric with partial polyethylene coating, is suitable as protective clothing against liquid chemicals and infectious agents. This durable isolation gown with knitted cuffs is classified as Personal Protective Equipment Category III and is liquid-and particle-repellent. It provides reliable protection against viruses, bacteria, and liquid aerosols.



**Tested to:**  
EN ISO 13688, EN 13034, EN 14325, EN 14126,  
EN ISO 6530, ISO 16603, ISO 16604, EN ISO 22610  
ISO/DIS 22611, EN ISO 22612, EN ISO 12947-2, EN  
863, EN ISO 13934-1, EN ISO 9073-4, EN ISO 13935-2

- PP spunbond with liquid-proof PE-part-coating (26 g/m<sup>2</sup>) on sleeves and front side
- Fastening tape at neck and tie belt at waist
- Knitted sleeve cuffs
- Latex-free
- Personal Protective Equipment Category III\*
- Type PB [6] according to EN 13034
- Type PB [6]-B according to EN 14126
- 10 pieces per bag

NEW  
Available  
Now



Item Number	Product	Available Size	Colour	Material & Grammage	Sales Unit (SU)
N11444	medi-inn CLASSIC Isolation Gown SELECT PPE Category III (with partial PE lamination)	139 x 139 cm	Yellow	PP Spunbond with PE part-coating 26 g/m <sup>2</sup>	10 x 10 pieces

\* According to Regulation (EU) 2016/425

## medi-inn CLASSIC Protective Gown with Full PE Coating

The medi-inn CLASSIC Protective Gown made of skin-friendly polypropylene non-woven fabric with full polyethylene coating is suitable for examinations and work in the medical field as well as hygiene prevention for nursing staff. The disposable gown with knitted cuffs can be used as a protective gown in basic care in the context of MRSA, MRGN and ESBL, noro and influenza infections, as well as for laboratory work and in the food industry.



Ties at the back of the neck and waist



Knitted cuffs on sleeves



Tested to:  
EN 14126:3002 / AC:2004, ISO 16603, ISO 16604,  
EN ISO 22610, ISO 22611, EN ISO 22612

- PP spunbond with full PE coating 32 g/m<sup>2</sup>
- Waist ties fasten at the sides
- Ties at the back of the neck and waist
- Knitted cuffs on sleeves
- Latex-free
- Category I personal protective equipment \*
- Class I medical device \*\*
- 10 pieces per bag



Item Number	Product	Available Size	Colour	Material & Grammage	Sales Unit (SU)
N11333-b	medi-inn CLASSIC Protective Gown (with full PE coating)	139 x 139 cm	Blue	PP spunbond with full PE coating 32 g/m <sup>2</sup>	10 x 10 Pieces
N11333-b-XL	medi-inn CLASSIC Protective Gown (with full PE coating)	140 x 160 cm	Blue	PP spunbond with full PE coating 32 g/m <sup>2</sup>	10 x 10 Pieces

\* According to Regulation (EU) 2016/425 \*\* According to Regulation (EU) 2017/745

## medi-inn CLASSIC Isolation Gown, PPE Cat. III

The medi-inn CLASSIC Isolation Gown made of skin-friendly polypropylene non-woven fabric with full polyethylene coating is suitable as protective clothing against liquid chemicals and infectious agents. The durable isolation gown with elastic cuffs is category III personal protective equipment and is liquid- and particle-repellent. It offers reliable protection against viruses, bacteria and liquid aerosols.



CE 0598  

- PP spunbond with full PE coating 42 g/m<sup>2</sup>
- Ultrasonically welded inner seams
- Waist tie belt fastens at the sides
- Fastening tape at neck and tie belt at waist
- Elastic sleeve cuffs
- Latex-free
- Personal Protective Equipment Category III \*
- Type PB [6] in accordance with EN 13034
- 10 pieces per bag



**Tested to:**  
EN ISO 13688, EN 13034, EN 14325, EN 14126,  
EN ISO 6530, ISO 16603, ISO 16604, EN ISO 22610  
ISO/DIS 22611, EN ISO 22612, EN ISO 12947-2, EN  
863, EN ISO 13934-1, EN ISO 9073-4, EN ISO 13935-2

Item Number	Product	Available Size	Colour	Material & Grammage	Sales Unit (SU)
N19062	medi-inn CLASSIC Isolation Gown, PPE Cat. III	139 x 139 cm	Orange	PP spunbond with full PE coating 42 g/m <sup>2</sup>	10 x 10 Pieces

\* According to Regulation (EU) 2016/425



[WWW.BODY-PRODUCTS.DE](http://WWW.BODY-PRODUCTS.DE)

Your contact:

**BODY PRODUCTS**  
Relax Pharma und Kosmetik GmbH

Alfred-Nobel-Straße 1-3  
50226 Frechen

Phone: +49 2234 91165-0  
Fax: +49 2234 91165-23  
Email: [info@body-products.de](mailto:info@body-products.de)