

medi-inn[®]
CLASSIC



Protective Gowns

Professionally Dressed.
Reliably Protected.

Professionally Dressed. Reliably Protected.

Hygiene awareness in the field of protective workwear has changed again in the recent past.

Our many years of experience in the production and marketing of medical protective gowns, the high quality of which is focused on the needs of our customers, means we are a reliable partner for workwear in hygiene-sensitive areas. Our protective gowns comply with European standards and regulations, and offer a variety of different styles and sizes for optimum freedom of movement.

As an expert in disposable hygiene products and protective clothing for hospitals, nursing and occupational safety, BODY PRODUCTS offers a carefully selected portfolio of visitor, protective and isolation gowns.



Hospital



Medical Practice



Dental Practice




Laboratory



Geriatric & Nursing Care

Special Features:



medi-inn® CLASSIC PP Visitor Gown



The medi-inn® CLASSIC PP Visitor Gown made of skin-friendly polypropylene non-woven fabric with elastic cuffs is breathable and is used as a disposable visitor gown for hygiene prevention for visitors and nursing staff or by the food industry.

- PP spunbond 16 g/m²
- Breathable
- Ties at the back of the neck and waist
- Elastic sleeve cuffs
- Latex-free
- Category I personal protective equipment *
- Class I medical device **
- 10 pieces wrapped in foil



Ties at the back of the neck and waist



Elastic sleeve cuffs

Item Number	Product	Available Size	Colour	Material & Grammage	Sales Unit (SU)
N10191-hb	medi-inn® CLASSIC PP Visitor Gown	115 x 137 cm	Light Blue	PP Spunbond 16 g/m ²	10 x 10 Pieces
N10191-w	medi-inn® CLASSIC PP Visitor Gown	115 x 137 cm	White	PP Spunbond 16 g/m ²	10 x 10 Pieces
N10191-hg	medi-inn® CLASSIC PP Visitor Gown	115 x 137 cm	Light Green	PP Spunbond 16 g/m ²	10 x 10 Pieces

* as per Regulation (EU) 2016/425 ** as per Regulation (EU) 2017/745

medi-inn® CLASSIC PP Visitor Gown



The medi-inn® CLASSIC PP Visitor Gown made of skin-friendly polypropylene non-woven fabric with knitted cuffs is hard-wearing and breathable. It is used as a disposable visitor gown for hygiene prevention for visitors and nursing staff or by the food industry.

- PP spunbond 24 g/m²
- Breathable
- Ties at the back of the neck and waist
- Knitted cuffs on sleeves
- Latex-free
- Category I personal protective equipment *
- Class I medical device **
- 10 pieces wrapped in foil



Item Number	Product	Available Size	Colour	Material & Grammage	Sales Unit (SU)
N11710-dg	medi-inn® CLASSIC PP Visitor Gown	136 x 140 cm	Dark Green	PP Spunbond 24 g/m ²	10 x 10 Pieces

medi-inn® CLASSIC Protective Gown with PE Part-Coating

The disposable gown with knitted cuffs can be used as a protective gown in basic care in the context of MRSA, MRGN and ESBL, noro and influenza infections, as well as for laboratory work and in the food industry.



Tested to:
EN 14126:2003/AC:2004
ISO 16603, ISO 16604,
EN ISO 22610, ISO 22611,
EN ISO 22612

- PP spunbond with liquid-tight PE part-coating on the sleeves and front
- Waist ties fasten at the sides
- Ties at the back of the neck and waist
- Knitted cuffs on sleeves
- Category I personal protective equipment *
- Class I medical device **
- Breathable and latex-free
- 10 pieces wrapped in foil



Item Number	Product	Available Size	Colour	Material & Grammage	Sales Unit (SU)
N11333-y	medi-inn® CLASSIC Protective Gown with PE part-coating	139 x 139 cm	Yellow	PP Spunbond with PE part-coating 31 g/m ²	10 x 10 Pieces
N11333-y-XL	medi-inn® CLASSIC Protective Gown with PE part-coating	140 x 160 cm	Yellow	PP Spunbond with PE part-coating 31 g/m ²	10 x 10 Pieces
N18344-y	medi-inn® CLASSIC Protective Gown with PE part-coating	136 x 136 cm	Yellow	PP Spunbond with PE part-coating 26 g/m ²	10 x 10 Pieces

* as per Regulation (EU) 2016/425 ** as per Regulation (EU) 2017/745

medi-inn® CLASSIC Protective Gown with Full PE Coating

The medi-inn® CLASSIC Protective Gown made of skin-friendly polypropylene non-woven fabric with full polyethylene coating is suitable for examinations and work in the medical field as well as hygiene prevention for nursing staff. The disposable gown with knitted cuffs can be used as a protective gown in basic care in the context of MRSA, MRGN and ESBL, noro and influenza infections, as well as for laboratory work and in the food industry.



Tested to:
EN 14126:3002 / AC:2004
ISO 16603, ISO 16604,
EN ISO 22610, ISO 22611,
EN ISO 22612

- PP spunbond with full PE coating 32 g/m²
- Waist ties fasten at the sides
- Ties at the back of the neck and waist
- Knitted cuffs on sleeves
- Latex-free
- Category I personal protective equipment *
- Class I medical device **
- 10 pieces wrapped in foil



Ties at the back
of the neck
and waist



Knitted cuffs
on sleeves

Item Number	Product	Available Size	Colour	Material & Grammage	Sales Unit (SU)
N11333-b	medi-inn® CLASSIC Protective Gown with full PE coating	139 x 139 cm	Blue	PP spunbond with full PE coating 32 g/m ²	10 x 10 Pieces
N11333-b-XL	medi-inn® CLASSIC Protective Gown with full PE coating	140 x 160 cm	Blue	PP spunbond with full PE coating 32 g/m ²	10 x 10 Pieces

* as per Regulation (EU) 2016/425 ** as per Regulation (EU) 2017/745

medi-inn® CLASSIC Isolation Gown, PPE Cat. III

The medi-inn® CLASSIC Isolation Gown made of skin-friendly polypropylene non-woven fabric with full polyethylene coating is suitable as protective clothing against liquid chemicals and infectious agents. The durable isolation gown with elastic cuffs is category III personal protective equipment and is liquid- and particle-repellent. It offers reliable protection against viruses, bacteria and liquid aerosols.

CE 0598  

Tested to:
EN ISO 13688, EN 13034,
EN 14325, EN 14126,
EN ISO 6530, ISO 16603,
ISO 16604, EN ISO 22610
ISO/DIS 22611, EN ISO 22612,
EN ISO 12947-2, EN 863,
EN ISO 13934-1, EN ISO 9073-4,
EN ISO 13935-2

- PP spunbond with full PE coating 42 g/m²
- Ultrasonically welded inner seams
- Waist tie belt fastens at the sides
- Fastening tape at neck and tie belt at waist
- Elastic sleeve cuffs
- Latex-free
- Category III personal protective equipment *
- Type PB [6] as per EN 13034
- 10 pieces wrapped in foil



Ties at the back
of neck
and waist



Elastic sleeve
cuffs



Waist tie
belt ties
at the side



View of outer seams

Ultrasonically
welded
inner seams



View of inner seams

Item Number	Product	Available Size	Colour	Material & Grammage	Sales Unit (SU)
N19062	medi-inn® CLASSIC Isolation Gown, PPE Cat. III	139 x 139 cm	Orange	PP spunbond with full PE coating 42 g/m ²	10 x 10 Pieces

* as per Regulation (EU) 2016/425



WWW.BODY-PRODUCTS.DE

Your contact:

BODY PRODUCTS

Relax Pharma und Kosmetik GmbH

Alfred-Nobel-Straße 1-3

50226 Frechen, Germany

Phone: +49 2234 91165-0

Fax: +49 2234 91165-23

Email: info@body-products.de