

medi-inn[®]
PRO

Protection

medi-inn[®] PRO Nitril CLINIC PERFECT

Nitrile Examination Gloves

For excellent protection against
viruses, bacteria, chemicals,
cytostatics and opioids.

 **BODY
PRODUCTS**

Developed for Professionals



Nitril CLINIC PERFECT Examination Gloves. Reliable Protection for Everyday Clinical Life.

**Powder-free, sensitive, tactile and pleasant to wear.
Highest functionality in use.**

Our examination gloves in the **medi-inn® PRO Nitril CLINIC PERFECT** product line have proved their worth in everyday clinical practice. Due to the special material composition, the durable examination gloves are suitable for challenging activities under extreme conditions. As well as being extremely pleasant to wear and offering good sensitivity, the focus is on protection against viruses, bacteria and chemicals as well as cytostatics.



Advantages in Clinical Use

- Highly tactile and sensitive grip.
- Pleasant wear comfort even during longer activities.
- Wide range of applications.
- Excellent protection against viruses, bacteria, chemicals and cytostatics.
- Suitable for users with a natural latex intolerance.



Advantages in Storage

- Optimised packaging dimensions ensure overall cost savings and an excellent price-performance ratio.
- An improved glove dispenser opening enables easy and quick removal of the gloves.

medi-inn® PRO Nitril CLINIC PERFECT WHITE

medi-inn® PRO Nitril CLINIC PERFECT BLUE

Our nitrile examination gloves have proven themselves in everyday clinical practice. Thanks to the excellent material composition, **medi-inn® PRO CLINIC PERFECT WHITE / BLUE** offer excellent protection against viruses, bacteria and chemicals as well as cytostatics despite their low wall thickness.



Product Properties	medi-inn® PRO Nitril CLINIC PERFECT WHITE	medi-inn® PRO Nitril CLINIC PERFECT BLUE
	Item No.: N12791-w	Item No.: N12791-b
	Available in white	Available in cobalt blue
	200 pieces per box 2,000 pieces per carton	
	Ambidextrous, powder-free (chlorinated), finger-textured	
	Length 240 mm	
	AQL 1.0	
	Complies with Medical Device Regulation (EU) 2017/745 (MD Class I) and PPE Regulation (EU) 2016/425 (PPE Category III)	
Wall Thickness (single wall)		
Fingers	0.07 mm	
Palm	0.06 mm	
Cuff	0.05 mm	



Protection against Chemicals, Cytostatics and an Opioid



18 Chemicals



15 Cytostatics / 1 Opioid



You will find detailed information in our technical data sheets

medi-inn® PRO Nitril CLINIC PERFECT SEMI LONG medi-inn® PRO Nitril CLINIC PERFECT LONG

Thanks to the greater wall thickness, our nitrile examination gloves in the **medi-inn® PRO CLINIC PERFECT SEMI LONG / LONG** line offer reinforced barrier properties. With extended cuffs that also provide additional protection for the forearm.



Product Properties	medi-inn® PRO Nitril CLINIC PERFECT SEMI LONG	medi-inn® PRO Nitril CLINIC PERFECT LONG
	Item No.: N13793	Item No.: N15091
	Available in white	Available in violet-blue
	150 pieces per box 1,500 pieces per carton	100 pieces per box 1,000 pieces per carton
	Ambidextrous, powder-free (chlorinated), finger-textured	
	Length 265 mm	Length 290 mm
	AQL 1.0	
	Complies with Medical Device Regulation (EU) 2017/745 (MD Class I) and PPE Regulation (EU) 2016/425 (PPE Category III)	
Wall Thickness (single wall)		
Fingers	0.12 mm	0.12 mm
Palm	0.07 mm	0.08 mm
Cuff	0.05 mm	0.06 mm



Protection against Chemicals and Cytostatics



16 Chemicals



12 Cytostatics

SEMI LONG



LONG



You will find detailed information in our technical data sheets

NEW

medi-inn[®]
PRO

medi-inn[®] PRO Nitril CLINIC PERFECT RECOVERY

The patented colloidal oatmeal inner coating of our nitrile examination glove provides a natural moisturising barrier between the glove and the skin. This can act as a preventative measure against irritated, dry and itchy skin. In particular, we recommend our **medi-inn[®] PRO Nitril CLINIC PERFECT RECOVERY** gloves to ensure that the hands of users who suffer from hypersensitivity of the skin due to frequent hand washing and the wearing of examination gloves stay soft and smooth.



U.S. Food and Drug Administration (FDA) approved skin protectant. Certified gluten-free glove that is also suitable for people with coeliac disease or who are sensitive to gluten.

Product Properties	medi-inn [®] PRO Nitril CLINIC PERFECT RECOVERY
	Item No.: N640528
	Available in grey-blue
	100 pieces per box / 1,000 pieces per carton
	Ambidextrous, powder-free (chlorinated), finger-textured
	COATS [®] inner coating
	Length 240 mm
	AQL 1.0
	Complies with Medical Device Regulation (EU) 2017/745 (MD Class I) and PPE Regulation (EU) 2016/425 (PPE Category III)
Wall Thickness (single wall)	
Fingers	0.09 mm
Palm	0.07 mm
Cuff	0.06 mm



Protection against Chemicals and Cytostatics



11 Chemicals



12 Cytostatics

You will find detailed information in our technical data sheets



Our Product Line medi-inn® PRO Nitril CLINIC PERFECT



Item Number	Product	Colour	Length	Contents of a Shipping Unit	Available Sizes
N12791-w	medi-inn® PRO Nitril CLINIC PERFECT WHITE	White	240 mm	10 x 200 pieces	XS - XL
N12791-b	medi-inn® PRO Nitril CLINIC PERFECT BLUE	Cobalt blue	240 mm	10 x 200 pieces	XS - XL
N13793	medi-inn® PRO Nitril CLINIC PERFECT SEMI LONG	White	265 mm	10 x 150 pieces	S - XL
N15091	medi-inn® PRO Nitril CLINIC PERFECT LONG	Violet-blue	290 mm	10 x 100 pieces	XS - XL
N640528	medi-inn® PRO Nitril CLINIC PERFECT RECOVERY	Grey-blue	240 mm	10 x 100 pieces	XS - XL