



Product Catalogue

For Hospitals, Care and Occupational Safety

WWW.BODY-PRODUCTS.DE

A WARM WELCOME

to the **BODY PRODUCTS** catalogue
for hospitals, care facilities and occupational safety.

Whether in hospitals, care facilities, medical practices, laboratories or in the field of occupational safety – hygienic working conditions and reliable protection are essential in your daily work. As experts in disposable hygiene products and personal protective equipment for hospitals, care and occupational safety, **BODY PRODUCTS** offers you practical, safe and convenient solutions for a wide range of requirements.

This catalogue provides you with a comprehensive overview of our product range – carefully selected and compiled in a user-friendly manner.

Our Brands

Systematic quality and variety

With our four brands, we offer you products for different applications and requirements:



medi-inn® CLASSIC

stands for proven standard products with tested quality at fair prices – ideal for everyday use in medical practices and care settings for healthcare professionals, patients, and caregivers.



medi-inn® PRO

is designed for healthcare professionals with higher demands on functionality and performance – developed for use in specialized medical and care environments.



medi-inn® GREEN

combines high product safety with more sustainable materials and environmentally-conscious thinking – for those who want to procure in a more ecologically responsible way.



WORK-INN®

offers high-quality products for occupational safety, catering and food processing – developed together with our partner PAPSTAR and designed according to HACCP standards for maximum safety in everyday working life.

Tested Quality

We attach great importance to the quality and performance of all products – which is why they are comprehensively tested. Our quality management system is certified according to DIN EN ISO 13485, and therefore meets the high requirements for medical devices.

The topic of sustainability is also firmly anchored in our corporate philosophy: we are certified according to the environmental standard DIN EN ISO 14001; our sustainability performance is regularly evaluated by EcoVadis, and we have documented our voluntary commitment in the German Sustainability Code (DNK).



WWW.BODY-PRODUCTS.DE

About BODY PRODUCTS

BODY PRODUCTS offers a carefully selected, comprehensive portfolio with more than 2,000 medical devices and healthcare items. Our portfolio covers the areas of protection, hygiene, diagnostics and incontinence care. This includes high-quality medical aids, ward and practice supplies as well as first aid equipment. Another key focus lies on disposable protective clothing. Innovative solutions for highly specialized medical applications – such as endoscopy, surgery and urology – complete our range.

The experts for medical devices, disposable hygiene products and personal protective equipment



Contents

Chapter 1 | Gloves

Pages 6–27

- Nitrile Examination and Protective Gloves
 - Nitrile Chemical Protection Gloves
 - Latex Examination and Protective Gloves
 - Vinyl Examination and Protective Gloves
 - Polyethylene Gloves
 - Cotton Tricot Gloves
 - Stainless Steel Dispenser for Disposable Gloves
 - Surgical Gloves
-

Chapter 2 | Protective Clothing

Pages 28–44

- Protective Gowns (Partially Coated, Fully Coated)
 - Visitor Gowns
 - Aprons
 - Apron Accessories
 - Medical Face Masks
 - FFP Respirators
 - Overshoes
 - Anti-Slip Socks
-

Chapter 3 | Ward Supplies

Pages 45–60

- Denture Containers
 - Patient Belongings Bags
 - Medication Administration
 - Medication Dispensers
 - Tablet Splitters
 - Urinal Bottles and Accessories
 - Bedpans
 - Kidney Dishes
 - Emesis Bags
 - Cosmetics
-

Chapter 4 | Diagnostic Aids

Pages 61–70

- Stethoscopes
 - Blood Pressure Monitors
 - Thermometers
 - Urine Specimen Container
 - Tourniquets
 - Telemetry Pouches
 - Wooden Tongue Depressors
 - Diagnostic Examination Lights
-

Chapter 5 | Nursing Care

Pages 71–81

- Waterless Body Washing
- Wash Gloves
- Hair Removal
- Oral Care
- Bibs
- Feeding Cups

Chapter 6 | Hygiene

Pages 82–86

- Examination Couch Rolls
 - Patient Transfer Sheets
 - Waste Disposal
-

Chapter 7 | First Aid

Pages 87–90

- Emergency Blankets
 - Instant Cold Packs
 - Reusable Hot / Cold Packs
-

Chapter 8 | Incontinence Care

Pages 91–100

- Pads
 - Shaped Pads
 - Incontinence Slips
 - Pants
 - Fixation Pants
 - Protective Underpads
 - Bed Protection
-

Chapter 9 | Urology

Pages 101–107

- Urine Bags
 - Urine Measurement System
 - Silicone Balloon Catheter
-

Chapter 10 | Surgery

Pages 108–110

- Sterile Specimen Retrieval Bags with Handle
 - Sterile Specimen Retrieval Bags without Handle
-

Chapter 11 | Operating Room Supplies

Pages 111–126

- Preoperative Preparation and Care
 - Patient Positioning and Fixation
 - Operating Table Covers
 - OP Floor Mats
 - Surgical Apparel
-

Chapter 12 | Dispenser Systems

Pages 127–136

- Wall-Mounted Dispenser
- Wall-Mounted Dispenser Accessories
- Stainless Steel Dispensers



01 Gloves

- Nitrile Examination and Protective Gloves
- Nitrile Chemical Protective Gloves
- Latex Examination and Protective Gloves
- Vinyl Examination and Protective Gloves
- Polyethylene Gloves
- Cotton Tricot Gloves
- Stainless Steel Dispensers for Disposable Gloves
- Surgical Gloves

medi-inn® CLASSIC Nitril BLUE SOFT

medi-inn® CLASSIC Nitril WHITE SOFT

Nitrile examination and protective gloves



- Ambidextrous
- Powder-free (chlorinated)
- Finger-textured
- Non-sterile
- Length: 240 mm
- Wall thickness: Fingers: 0.07 mm / Palm: 0.06 mm
- AQL 1.5
- Class I medical device in accordance with Regulation (EU) 2017/745
- Category III personal protective equipment in accordance with Regulation (EU) 2016/425
- Chemical protection class type B
- Suitable for food contact in accordance with Regulation (EC) 1935/2004



CE 2777



Item number	Size	Colour	Pieces per box	Boxes per carton
N11248-S	S	Blue	100 pieces	10 boxes
N11248-M	M	Blue	100 pieces	10 boxes
N11248-L	L	Blue	100 pieces	10 boxes
N11248-XL	XL	Blue	100 pieces	10 boxes
N11252-S	S	White	100 pieces	10 boxes
N11252-M	M	White	100 pieces	10 boxes
N11252-L	L	White	100 pieces	10 boxes
N11252-XL	XL	White	100 pieces	10 boxes

You will find further information on the standards and recommendations applied to chemicals and cytostatics in our data sheets at www.body-products.de

medi-inn® PRO Nitril CLINIC PERFECT BLUE

medi-inn® PRO Nitril CLINIC PERFECT WHITE



Nitrile examination and protective gloves

- Ambidextrous
- Powder-free (chlorinated)
- Finger-textured
- Non-sterile
- Length: 240 mm
- Wall thickness: Fingers: 0.07 mm / Palm: 0.06 mm
- AQL 1.0
- With patented LOW DERMA™ technology
- Class I medical device in accordance with Regulation (EU) 2017/745
- Category III personal protective equipment in accordance with Regulation (EU) 2016/425
- Chemical protection class type B
- Suitable for food contact in accordance with Regulation (EC) 1935/2004
- Protection against chemicals, cytostatics and an opioid*



Item number	Size	Colour	Pieces per box	Boxes per carton
N12791-b-XS	XS	Cobalt blue	200 pieces	10 boxes
N12791-b-S	S	Cobalt blue	200 pieces	10 boxes
N12791-b-M	M	Cobalt blue	200 pieces	10 boxes
N12791-b-L	L	Cobalt blue	200 pieces	10 boxes
N12791-b-XL	XL	Cobalt blue	200 pieces	10 boxes
N12791-w-XS	XS	White	200 pieces	10 boxes
N12791-w-S	S	White	200 pieces	10 boxes
N12791-w-M	M	White	200 pieces	10 boxes
N12791-w-L	L	White	200 pieces	10 boxes
N12791-w-XL	XL	White	200 pieces	10 boxes

*You will find further information on the standards and recommendations applied to chemicals and cytostatics in our data sheets at www.body-products.de

medi-inn® CLASSIC Nitril BLUE

Nitrile examination and protective gloves



- Ambidextrous
- Powder-free (chlorinated)
- Finger-textured
- Non-sterile
- Length: 240 mm
- Wall thickness: Fingers: 0.08 mm / Palm: 0.06 mm
- AQL 1.5
- Class I medical device in accordance with Regulation (EU) 2017/745
- Category III personal protective equipment in accordance with Regulation (EU) 2016/425
- Chemical protection class type B
- Suitable for food contact in accordance with Regulation (EC) 1935/2004



CE 2777



Item number	Size	Colour	Pieces per box	Boxes per carton
N10419-S	S	Blue	200 pieces	10 boxes
N10419-M	M	Blue	200 pieces	10 boxes
N10419-L	L	Blue	200 pieces	10 boxes
N10419-XL	XL	Blue	200 pieces	10 boxes

You will find further information on the standards and recommendations applied to chemicals and cytostatics in our data sheets at www.body-products.de

medi-inn® CLASSIC Nitril BLUE PLUS

medi-inn® CLASSIC Nitril WHITE PLUS

Nitrile examination and protective gloves



- Ambidextrous
- Powder-free (chlorinated)
- Finger-textured
- Non-sterile
- Length: 240 mm
- Wall thickness: Fingers: 0.10 mm / Palm: 0.07 mm
- AQL 1.5
- Class I medical device in accordance with Regulation (EU) 2017/745
- Category III personal protective equipment in accordance with Regulation (EU) 2016/425
- Chemical protection class type B
- Suitable for food contact in accordance with Regulation (EC) 1935/2004



CE 2777



Item number	Size	Colour	Pieces per box	Boxes per carton
N10111-XS	XS (on enquiry)	Blue	100 pieces	10 boxes
N10111-S	S	Blue	100 pieces	10 boxes
N10111-M	M	Blue	100 pieces	10 boxes
N10111-L	L	Blue	100 pieces	10 boxes
N10111-XL	XL	Blue	100 pieces	10 boxes
N10111-2X	XXL	Blue	100 pieces	10 boxes
N17322-S	S	White	100 pieces	10 boxes
N17322-M	M	White	100 pieces	10 boxes
N17322-L	L	White	100 pieces	10 boxes
N17322-XL	XL	White	100 pieces	10 boxes

You will find further information on the standards and recommendations applied to chemicals and cytostatics in our data sheets at www.body-products.de

medi-inn® PRO Nitril CLINIC PERFECT RECOVERY



Nitrile examination and protective gloves

- Ambidextrous
- Powder-free (chlorinated)
- Finger-textured
- COATS® Inner Coating
- Non-sterile
- Length 240 mm
- Wall thickness: 0.09 mm / Palm: 0.07 mm
- AQL 1.0
- Dermatest® tested
- Class I medical device in accordance with Regulation (EU) 2017/745
- Category III personal protective equipment in accordance with Regulation (EU) 2016/425
- Suitable for short-term food contact in accordance with Regulation (EC) 1935/2004



Item number	Size	Colour	Pieces per box	Boxes per carton
N640528-XS	XS	Grey-Blue	100 pieces	10 boxes
N640528-S	S	Grey-Blue	100 pieces	10 boxes
N640528-M	M	Grey-Blue	100 pieces	10 boxes
N640528-L	L	Grey-Blue	100 pieces	10 boxes
N640528-XL	XL	Grey-Blue	100 pieces	10 boxes

You will find further information on the standards and recommendations applied to chemicals and cytostatics in our data sheets at www.body-products.de

medi-inn® CLASSIC Nitril



Nitrile examination and protective gloves

- Available in 10 different colours
- Ambidextrous
- Powder-free (chlorinated)
- Finger-textured
- Non-sterile
- Length: 240 mm
- Wall thickness: Fingers: 0.09 - 0.10 mm / Palm: 0.07 mm
- AQL 1.5
- Class I medical device in accordance with Regulation (EU) 2017/745
- Category III personal protective equipment in accordance with Regulation (EU) 2016/425
- Chemical protection class type B
- Suitable for food contact in accordance with Regulation (EC) 1935/2004



Available Colours



Item number	Product name	Size	Colour	Pieces per box	Boxes per carton
N11521-XS	medi-inn® CLASSIC Nitril BLACK PLUS	XS	Black	100 pieces	10 boxes
N11521-S	medi-inn® CLASSIC Nitril BLACK PLUS	S	Black	100 pieces	10 boxes
N11521-M	medi-inn® CLASSIC Nitril BLACK PLUS	M	Black	100 pieces	10 boxes
N11521-L	medi-inn® CLASSIC Nitril BLACK PLUS	L	Black	100 pieces	10 boxes
N11521-XL	medi-inn® CLASSIC Nitril BLACK PLUS	XL	Black	100 pieces	10 boxes
N15230-XS	medi-inn® CLASSIC Nitril GREEN PLUS	XS	Green	100 pieces	10 boxes
N15230-S	medi-inn® CLASSIC Nitril GREEN PLUS	S	Green	100 pieces	10 boxes
N15230-M	medi-inn® CLASSIC Nitril GREEN PLUS	M	Green	100 pieces	10 boxes
N15230-L	medi-inn® CLASSIC Nitril GREEN PLUS	L	Green	100 pieces	10 boxes
N15230-XL	medi-inn® CLASSIC Nitril GREEN PLUS	XL	Green	100 pieces	10 boxes
N18747-XS	medi-inn® CLASSIC Nitril MINT	XS	Mint	100 pieces	10 boxes
N18747-S	medi-inn® CLASSIC Nitril MINT	S	Mint	100 pieces	10 boxes
N18747-M	medi-inn® CLASSIC Nitril MINT	M	Mint	100 pieces	10 boxes
N18747-L	medi-inn® CLASSIC Nitril MINT	L	Mint	100 pieces	10 boxes
N18747-XL	medi-inn® CLASSIC Nitril MINT	XL	Mint	100 pieces	10 boxes

Item number	Product name	Size	Colour	Pieces per box	Boxes per carton
N18749-XS	medi-inn® CLASSIC Nitril MAGENTA	XS	Magenta	100 pieces	10 boxes
N18749-S	medi-inn® CLASSIC Nitril MAGENTA	S	Magenta	100 pieces	10 boxes
N18749-M	medi-inn® CLASSIC Nitril MAGENTA	M	Magenta	100 pieces	10 boxes
N18749-L	medi-inn® CLASSIC Nitril MAGENTA	L	Magenta	100 pieces	10 boxes
N18749-XL	medi-inn® CLASSIC Nitril MAGENTA	XL	Magenta	100 pieces	10 boxes
N18750-XS	medi-inn® CLASSIC Nitril PINK PLUS	XS	Pink	100 pieces	10 boxes
N18750-S	medi-inn® CLASSIC Nitril PINK PLUS	S	Pink	100 pieces	10 boxes
N18750-M	medi-inn® CLASSIC Nitril PINK PLUS	M	Pink	100 pieces	10 boxes
N18750-L	medi-inn® CLASSIC Nitril PINK PLUS	L	Pink	100 pieces	10 boxes
N18750-XL	medi-inn® CLASSIC Nitril PINK PLUS	XL	Pink	100 pieces	10 boxes
N18751-XS	medi-inn® CLASSIC Nitril PURPLE	XS	Purple	100 pieces	10 boxes
N18751-S	medi-inn® CLASSIC Nitril PURPLE	S	Purple	100 pieces	10 boxes
N18751-M	medi-inn® CLASSIC Nitril PURPLE	M	Purple	100 pieces	10 boxes
N18751-L	medi-inn® CLASSIC Nitril PURPLE	L	Purple	100 pieces	10 boxes
N18751-XL	medi-inn® CLASSIC Nitril PURPLE	XL	Purple	100 pieces	10 boxes
N18752-XS	medi-inn® CLASSIC Nitril RED PLUS	XS	Red	100 pieces	10 boxes
N18752-S	medi-inn® CLASSIC Nitril RED PLUS	S	Red	100 pieces	10 boxes
N18752-M	medi-inn® CLASSIC Nitril RED PLUS	M	Red	100 pieces	10 boxes
N18752-L	medi-inn® CLASSIC Nitril RED PLUS	L	Red	100 pieces	10 boxes
N18752-XL	medi-inn® CLASSIC Nitril RED PLUS	XL	Red	100 pieces	10 boxes
N18753-XS	medi-inn® CLASSIC Nitril YELLOW	XS	Yellow	100 pieces	10 boxes
N18753-S	medi-inn® CLASSIC Nitril YELLOW	S	Yellow	100 pieces	10 boxes
N18753-M	medi-inn® CLASSIC Nitril YELLOW	M	Yellow	100 pieces	10 boxes
N18753-L	medi-inn® CLASSIC Nitril YELLOW	L	Yellow	100 pieces	10 boxes
N18753-XL	medi-inn® CLASSIC Nitril YELLOW	XL	Yellow	100 pieces	10 boxes
N18755-XS	medi-inn® CLASSIC Nitril ORANGE	XS	Orange	100 pieces	10 boxes
N18755-S	medi-inn® CLASSIC Nitril ORANGE	S	Orange	100 pieces	10 boxes
N18755-M	medi-inn® CLASSIC Nitril ORANGE	M	Orange	100 pieces	10 boxes
N18755-L	medi-inn® CLASSIC Nitril ORANGE	L	Orange	100 pieces	10 boxes
N18755-XL	medi-inn® CLASSIC Nitril ORANGE	XL	Orange	100 pieces	10 boxes
N18886-XS	medi-inn® CLASSIC Nitril APPLE GREEN	XS	Apple green	100 pieces	10 boxes
N18886-S	medi-inn® CLASSIC Nitril APPLE GREEN	S	Apple green	100 pieces	10 boxes
N18886-M	medi-inn® CLASSIC Nitril APPLE GREEN	M	Apple green	100 pieces	10 boxes
N18886-L	medi-inn® CLASSIC Nitril APPLE GREEN	L	Apple green	100 pieces	10 boxes
N18886-XL	medi-inn® CLASSIC Nitril APPLE GREEN	XL	Apple green	100 pieces	10 boxes

You will find further information on the standards and recommendations applied to chemicals and cytostatics in our data sheets at www.body-products.de



medi-inn® GREEN Nitril

Biodegradable nitrile examination and protective gloves

- Biodegradable
- Ambidextrous
- Powder-free (chlorinated)
- Finger-textured
- Non-sterile
- Length: 240 mm
- Wall thickness: Fingers: 0.10 mm / Palm: 0.07 mm
- AQL 1.5
- Class I medical device in accordance with Regulation (EU) 2017/745
- Category III personal protective equipment in accordance with Regulation (EU) 2016/425
- Suitable for food contact in accordance with Regulation (EC) 1935/2004
- The biodegradability was tested to standards ASTM D5526 and ASTM D5511



Item number	Size	Colour	Pieces per box	Boxes per carton
N19034-XS	XS	Green	100 pieces	10 boxes
N19034-S	S	Green	100 pieces	10 boxes
N19034-M	M	Green	100 pieces	10 boxes
N19034-L	L	Green	100 pieces	10 boxes
N19034-XL	XL	Green	100 pieces	10 boxes

You will find further information on the standards and recommendations applied to chemicals and cytostatics in our data sheets at www.body-products.de

WORK-INN® Nitril Blue Extra

Nitrile examination and protective gloves



- Ambidextrous
- Powder-free (chlorinated)
- Textured
- Non-sterile
- Length: 240 mm
- Wall thickness: Fingers: 0.11 mm / Palm: 0.09 mm
- AQL 1.5
- Class I medical device in accordance with Regulation (EU) 2017/745
- Category III personal protective equipment in accordance with Regulation (EU) 2016/425
- Chemical protection class type B
- Suitable for food contact in accordance with Regulation (EC) 1935/2004



Item number	Size	Colour	Pieces per box	Boxes per carton
N18894-S	S	Blue	100 pieces	10 boxes
N18894-M	M	Blue	100 pieces	10 boxes
N18894-L	L	Blue	100 pieces	10 boxes
N18894-XL	XL	Blue	100 pieces	10 boxes
N18894-XXL	XXL	Blue	100 pieces	10 boxes

You will find further information on the standards and recommendations applied to chemicals and cytostatics in our data sheets at www.body-products.de

medi-inn® CLASSIC Nitril BLUE LONG

Nitrile examination and protective gloves



- Ambidextrous
- Powder-free (chlorinated)
- Textured
- Non-sterile
- Length: 300 mm
- Wall thickness: Fingers: 0.14 mm / Palm: 0.10 mm
- AQL 1.5
- Class I medical device in accordance with Regulation (EU) 2017/745
- Category III personal protective equipment in accordance with Regulation (EU) 2016/425
- Chemical protection class type B
- Suitable for food contact in accordance with Regulation (EC) 1935/2004



Item number	Size	Colour	Pieces per box	Boxes per carton
N11872-S	S	Blue	100 pieces	10 boxes
N11872-M	M	Blue	100 pieces	10 boxes
N11872-L	L	Blue	100 pieces	10 boxes
N11872-XL	XL	Blue	100 pieces	10 boxes

You will find further information on the standards and recommendations applied to chemicals and cytostatics in our data sheets at www.body-products.de

WORK-INN® Nitril High Risk Blue

Nitrile chemical protective gloves



- Ambidextrous
- Powder-free (chlorinated)
- Textured
- Non-sterile
- Length: 300 mm
- Wall thickness: Fingers: 0.32 mm / Palm: 0.26 mm
- AQL 1.5
- Class I medical device in accordance with Regulation (EU) 2017/745
- Category III personal protective equipment in accordance with Regulation (EU) 2016/425
- Chemical protection class type A
- Suitable for food contact in accordance with Regulation (EC) 1935/2004



Item number	Size	Colour	Pieces per box	Boxes per carton
N18897-S	S	Blue	50 pieces	10 boxes
N18897-M	M	Blue	50 pieces	10 boxes
N18897-L	L	Blue	50 pieces	10 boxes
N18897-XL	XL	Blue	50 pieces	10 boxes
N18897-XXL	XXL	Blue	50 pieces	10 boxes

You will find further information on the standards and recommendations applied to chemicals and cytostatics in our data sheets at www.body-products.de

medi-inn® CLASSIC Latex WHITE GRIP

Latex examination and protective gloves



- Ambidextrous
- Powder free (polymer coated)
- Textured
- Non-sterile
- Length: 240 mm
- Wall thickness: Fingers: 0.12 mm / palm: 0.11 mm
- AQL 1.5
- Class I medical device in accordance with Regulation (EU) 2017/745
- Category III personal protective equipment in accordance with Regulation (EU) 2016/425
- Chemical protection class type B
- Suitable for food contact in accordance with Regulation (EC) 1935/2004



Item number	Size	Colour	Pieces per box	Boxes per carton
N10103-XS	XS (on enquiry)	White (natural)	100 pieces	10 boxes
N10103-S	S	White (natural)	100 pieces	10 boxes
N10103-M	M	White (natural)	100 pieces	10 boxes
N10103-L	L	White (natural)	100 pieces	10 boxes
N10103-XL	XL	White (natural)	100 pieces	10 boxes

You will find further information on the standards and recommendations applied to chemicals and cytostatics in our data sheets at www.body-products.de

WORK-INN® Latex Comfort



Latex examination and protective gloves

- Ambidextrous
- Powder-free (chlorinated)
- Textured
- Non-sterile
- Length: 240 mm
- Wall thickness: Fingers: 0.12 mm / Palm: 0.11 mm
- AQL 1.5
- Class I medical device in accordance with Regulation (EU) 2017/745
- Category III personal protective equipment in accordance with Regulation (EU) 2016/425
- Chemical protection class type B
- Suitable for food contact in accordance with Regulation (EC) 1935/2004



Item number	Size	Colour	Pieces per box	Boxes per carton
N18899-S	S	White (natural)	100 pieces	10 boxes
N18899-M	M	White (natural)	100 pieces	10 boxes
N18899-L	L	White (natural)	100 pieces	10 boxes
N18899-XL	XL	White (natural)	100 pieces	10 boxes

You will find further information on the standards and recommendations applied to chemicals and cytostatics in our data sheets at www.body-products.de



medi-inn® CLASSIC Latex PREMIUM

Latex examination and protective gloves

- Ambidextrous
- Powder-free (chlorinated)
- Textured
- Non-sterile
- Length: 240 mm
- Wall thickness: Fingers: 0.15 mm / Palm: 0.13 mm
- AQL 1.5
- Class I medical device in accordance with Regulation (EU) 2017/745
- Category III personal protective equipment in accordance with Regulation (EU) 2016/425
- Chemical protection class type B
- Suitable for food contact in accordance with Regulation (EC) 1935/2004



Item number	Size	Colour	Pieces per box	Boxes per carton
100376-S	S	White (natural)	100 pieces	10 boxes
100376-M	M	White (natural)	100 pieces	10 boxes
100376-L	L	White (natural)	100 pieces	10 boxes
100376-XL	XL	White (natural)	100 pieces	10 boxes

You will find further information on the standards and recommendations applied to chemicals and cytostatics in our data sheets at www.body-products.de

medi-inn® CLASSIC Vinyl COMFORT



Vinyl examination and protective gloves

- Ambidextrous
- Powder free (polymer coated)
- Smooth
- Non-sterile
- Length: 240 mm
- Wall thickness: Fingers: 0.08 mm / Palm: 0.07 mm
- AQL 1.5
- Class I medical device in accordance with Regulation (EU) 2017/745
- Category I personal protective equipment in accordance with Regulation (EU) 2016/425
- Suitable for short-term contact with non-greasy foods in accordance with Regulation (EC) 1935/2004



Item number	Size	Colour	Pieces per box	Boxes per carton
N10117-S	S	Transparent	100 pieces	10 boxes
N10117-M	M	Transparent	100 pieces	10 boxes
N10117-L	L	Transparent	100 pieces	10 boxes
N10117-XL	XL	Transparent	100 pieces	10 boxes

You will find further information on the standards and recommendations applied to chemicals and cytostatics in our data sheets at www.body-products.de

medi-inn® CLASSIC Vinyl STRETCH

Vinyl examination and protective gloves



- Ambidextrous
- Powder free (polymer coated)
- Smooth
- Non-sterile
- Length: 240 mm
- Wall thickness: Fingers: 0.08 mm / palm: 0.07 mm
- AQL 1.5
- Class I medical device in accordance with Regulation (EU) 2017/745
- Category I personal protective equipment in accordance with Regulation (EU) 2016/425
- Suitable for short-term contact with non-greasy foods in accordance with Regulation (EC) 1935/2004



Item number	Size	Colour	Pieces per box	Boxes per carton
N10125-S	S	Transparent	100 pieces	10 boxes
N10125-M	M	Transparent	100 pieces	10 boxes
N10125-L	L	Transparent	100 pieces	10 boxes
N10125-XL	XL	Transparent	100 pieces	10 boxes

You will find further information on the standards and recommendations applied to chemicals and cytostatics in our data sheets at www.body-products.de

PE Gloves

Protective gloves made of polyethylene

- Ambidextrous
- Powder free
- Hammered surface
- Non-sterile
- Wall thickness 0.012 mm
- Category I personal protective equipment in accordance with Regulation (EU) 2016/425
- Suitable for short-term food contact in accordance with Regulation (EC) 1935/2004
- Packed in a polybag of 100 pieces




Item number	Size	Dimensions	Colour	Pieces per bag	Bags per carton
N11110-STD	M	Length 280 mm, width 230 mm	Transparent	100 pieces	100 bags
N11112-STD	L	Length 300 mm, width 240 mm	Transparent	100 pieces	100 bags

medi-inn® CLASSIC Cotton Tricot Gloves



Custom-fit undergloves, protective gloves and allergy gloves made of 100% cotton

- 100% cotton
- Lined
- Breathable
- Non-sterile
- With sidewall and inserted thumb
- Very good sensitivity, perfect fit thanks to individual sizes
- Category I personal protective equipment in accordance with Regulation (EU) 2016/425
- Packed in foil bags of 12 pairs



Item number	Size	Colour	Pairs per bag	Bags per carton
N640622-7	7 / S	Pure white	12 pairs	50 bags
N640622-8	8 / M	Pure white	12 pairs	50 bags
N640622-9	9 / L	Pure white	12 pairs	50 bags
N640622-10	10 / XL	Pure white	12 pairs	50 bags
N640622-11	11 / XXL	Pure white	12 pairs	50 bags
N640622-12	12 / XXXL	Pure white	12 pairs	50 bags



medi-inn® CLASSIC Stainless Steel Dispensers for Disposable Gloves

Stainless steel dispensers for the practical storage of glove boxes and the quick removal of disposable gloves

- Dispenser for storing one glove box
- For wall mounting
- Made of stainless steel
- Including fittings
- Colour: Steel grey
- Suitable for standard glove boxes
Size XL suitable for glove box heights between 120 and 135 mm



Item number	Product name	Dimensions stainless steel dispenser Width / height / depth	Storage of	SU	Pieces per carton	Cartons per pallet
N12797-XL	Stainless steel dispenser XL	252 x 135 x 80 mm	1 glove box	1 piece	9 pieces	24

medi-inn® CLASSIC Multi-Box Stainless Steel Dispensers for Disposable Gloves



Stainless steel dispenser for the practical storage of glove boxes and the quick removal of disposable gloves

- Multi-Box dispenser for storing up to three or four glove boxes
- For wall mounting
- Made of stainless steel
- Including fittings
- Colour: Steel grey
- Suitable for standard glove boxes

Item no. N13218 suitable for glove box widths between 240 and 255 mm
Item no. N14097 suitable for glove box widths between 190 and 220 mm
Item no. N13219 suitable for glove box widths between 240 and 255 mm



Item example N13219

Item number	Product name	Dimensions stainless steel dispenser Width / height / depth	Storage of	SU	Pieces per carton	Cartons per pallet
N13218	3-way stainless steel dispenser	255 x 390 x 80 mm	3 glove boxes	1 piece	9 pieces	8
N14097	3-way stainless steel dispenser, narrow	220 x 350 x 70 mm	3 glove boxes	1 piece	9 pieces	8
N13219	4-way stainless steel dispenser	255 x 515 x 80 mm	4 glove boxes	1 piece	9 pieces	6



ulma® Surgical Gloves

Full protection for different requirements

Latex Surgical Gloves

Excellent biological barrier protection. Latex softness provides superior comfort, dexterity & flexibility.

Powdered | Powder free | Polymer coated

Chloroprene Surgical Gloves

Super softness & elasticity without risk of Type I allergy.

Powder free | Polymer coated

Polyisoprene Surgical Gloves

Top performing synthetic surgical gloves with super softness & elasticity without risk of Type I allergy.

Powder free | Polymer coated

Flexylon Surgical Gloves

Top quality without microchannels, featuring the highest AQL of 0.1. The accelerator-free latex formula prevents both Type I and Type IV allergies.

Powder free | Polymer coated

ulma® Product Range

Latex powdered

- ULMA® OP – General surgery

Latex powder free

- MAXITEX® PF – General surgery
- SENSIFLEX® PLUS – General surgery
- MAXITEX® UG – Undergloves
- MAXITEX® DUPLEX PF - Higher risk or intensive surgery
- MAXITEX® NEURO – Microsurgery
- MAXITEX® ORTHO – Orthopaedics

Chloroprene

- NUZONE® X2E – General surgery

Polyisoprene

- FUSIONE® – General surgery

Flexylon

- FINESSIS ZERO® – General surgery
- FINESSIS CORIUM® – Microsurgery and undergloves
- FINESSIS AEGIS® – General, orthopaedic, trauma and high-risk



Our Partner for Surgical Gloves



Ulma International GmbH is the European sales hub for the world's largest surgical glove manufacturer. Leveraging advanced technology and sustainable production practices, Ulma produces over 160 million surgical gloves every month. A strong commitment to quality, safety, and sustainability shapes all operations – from water-saving processes and green energy use to 92% recycled packaging. Ulma also takes responsibility beyond production, with a clear focus on inclusivity, community engagement, and workplace safety. Supported by a comprehensive product portfolio and outstanding customer service, Ulma International GmbH actively contributes to advancing healthcare standards in Europe and worldwide.



**CLICK HERE FOR THE COMPLETE
PRODUCT BROCHURE.**



Ulma International GmbH

Alfred-Nobel-Str. 5 | 50226 Frechen | Germany

☎ +49 731 926043-0

✉ info@ulma.de

www.ulma.de



Protective Clothing

- Protective Gowns (Partially Coated, Fully Coated)
- Visitor Gowns
- Aprons
- Apron Accessories
- Medical Face Masks
- FFP Respirators
- Overshoes
- Anti-Slip Socks



medi-inn® CLASSIC Protective Gown with Partial PE Coating

Partially coated protective gown

- PP spunbond with liquid-proof partial PE coating on the sleeves and front
- Knitted cuffs on sleeves
- Breathable
- Waist ties fasten at the sides
- Ties at the back of the neck and waist
- Latex-free
- Non-sterile
- Class I medical device in accordance with Regulation (EU) 2017/745
- Category I personal protective equipment in accordance with Regulation (EU) 2016/425



Item number	Material	Colour	Basis weight	Length x width	Multipack	Pieces per carton	Cartons per pallet
N11333-y	PP spunbond with liquid-proof partial PE Coating	Yellow	31 g / m ²	139 x 139 cm	10 pieces	100 pieces	18
N11333-y-XL	PP spunbond with liquid-proof partial PE Coating	Yellow	31 g / m ²	140 x 160 cm	10 pieces	100 pieces	12
N18344-y	PP spunbond with liquid-proof partial PE Coating	Yellow	26 g / m ²	136 x 136 cm	10 pieces	100 pieces	24

Tested in accordance with EN 14126:2003/AC:2004, ISO 16603, ISO 16604, EN ISO 22610, ISO 22611, EN ISO 22612



medi-inn® CLASSIC Protective Gown with Full PE Coating

Fully coated protective gown

- PP spunbond with full PE coating 32 g / m²
- Knitted cuffs on sleeves
- Liquid-proof
- Waist ties fasten at the sides
- Ties at the back of the neck and waist
- Latex free
- Non-sterile
- Class I medical device in accordance with Regulation (EU) 2017/745
- Category I personal protective equipment in accordance with Regulation (EU) 2016/425



Item number	Material	Colour	Basis weight	Length x width	Multipack	Pieces per carton	Cartons per pallet
N11333-b	PP spunbond with liquid-proof full PE Coating	Blue	32 g / m ²	139 x 139 cm	10 pieces	100 pieces	16
N11333-b-XL	PP spunbond with liquid-proof full PE Coating	Blue	32 g / m ²	140 x 160 cm	10 pieces	100 pieces	9

Tested in accordance with EN 14126:2003/AC:2004, ISO 16603, ISO 16604, EN ISO 22610, ISO 22611, EN ISO 22612



medi-inn® CLASSIC Isolation Gown PPE Cat. III

Fully coated protective gown

- PP spunbond with full PE coating 42 g / m²
- Elastic sleeve cuffs
- Liquid-proof
- Ultrasonically welded inner seams
- Waist tie belt fastens at the sides
- Fastening tape at neck and tie belt at waist
- Latex free
- Non-sterile
- Type PB [6] in accordance with EN 13034
- Category III personal protective equipment in accordance with Regulation (EU) 2016/425



CE 0598



Item number	Material	Colour	Basis weight	Length x width	Multipack	Pieces per carton	Cartons per pallet
N19062	PP spunbond with liquid-proof full PE Coating	Orange	42 g / m ²	139 x 139 cm	10 pieces	100 pieces	12

Tested to: EN ISO 13688, EN 13034, EN 14325, EN 14126, EN ISO 6530, ISO 16603, ISO 16604, EN ISO 22610, ISO/DIS 22611, EN ISO 22612, EN ISO 12947-2, EN 863, EN ISO 13934-1, EN ISO 9073-4, EN ISO13935-2



medi-inn® CLASSIC PP Visitor Gown

Non-coated protective gown

- PP spunbond 24 g / m²
- Knitted cuffs on sleeves
- Breathable
- Ties at the back of the neck and waist
- Latex-free
- Non-sterile
- Class I medical device in accordance with Regulation (EU) 2017/745
- Category I personal protective equipment in accordance with Regulation (EU) 2016/425



Item number	Material	Colour	Basis weight	Length x width	Multipack	Pieces per carton	Cartons per pallet
N11710-dg	PP spunbond	Dark green	24 g / m ²	136 x 140 cm	10 pieces	100 pieces	15

medi-inn® CLASSIC PE Apron, PLUS

Non-sterile PE apron, for single use only



- Made of polyethylene (PE)
- Hammered
- Liquid-proof
- Neck cut-out and waistbands
- Latex-free
- Non-sterile
- Category I personal protective equipment in accordance with Regulation (EU) 2016/425
- Suitable for food contact in accordance with Regulation (EC) 1935/2004. Tested in accordance with LFGB and Regulation (EU) 10/2011 for short-term food contact



Item number	Material	Colour	Material thickness	Length x width	Design	Multipack	Pieces per carton	Cartons per pallet
703040-b	Polyethylene (PE)	Blue	17 µm	125 x 75 cm	In bags with dispenser opening	100 pieces	1000 pieces	24



medi-inn® CLASSIC PE Apron, PLUS

Non-sterile PE apron, for single use only

- Made of polyethylene (PE)
- Hammered
- Liquid-proof
- Neck cut-out and waistbands
- Latex-free
- Non-sterile
- Category I personal protective equipment in accordance with Regulation (EU) 2016/425
- Suitable for food contact in accordance with Regulation (EC) 1935/2004. Tested in accordance with LFGB and Regulation (EU) 10/2011 for short-term food contact



Item number	Material	Colour	Material thickness	Length x width	Design	Multipack	Pieces per carton	Cartons per pallet
703040-w	Polyethylene (PE)	White	17 µm	125 x 75 cm	In bags with dispenser opening	100 pieces	1000 pieces	32
703041-w	Polyethylene (PE)	White	17 µm	140 x 75 cm	In bags with dispenser opening	100 pieces	1000 pieces	32
703042-w	Polyethylene (PE)	White	17 µm	160 x 75 cm	In bags with dispenser opening	100 pieces	1000 pieces	24

medi-inn® CLASSIC PE Apron, EXTRA STRONG

Non-sterile PE apron, for single use only



- Made of polyethylene (PE)
- Smooth
- Liquid-proof
- Neck cut-out and waistbands
- Non-sterile
- Category I personal protective equipment in accordance with Regulation (EU) 2016/425
- Suitable for food contact in accordance with Regulation (EG) 1935/2004. Tested in accordance with LFGB and Regulation (EU) 10/2011 for short-term food contact



Item number	Material	Colour	Material thickness	Length x width	Design	Multipack	Pieces per carton	Cartons per pallet
703043-b	Polyethylene (PE)	Blue	47 µm	120 x 75 cm	Blocked	50 pieces	500 pieces	24

medi-inn® CLASSIC PE Apron, STANDARD

Non-sterile PE aprons, for single use only



- Made of polyethylene (PE)
- Hammered
- Liquid-proof
- Neck cut-out and waistbands
- Latex-free
- Non-sterile
- Category I personal protective equipment in accordance with Regulation (EU) 2016/425
- Suitable for food contact in accordance with Regulation (EG) 1935/2004. Tested in accordance with LFGB and Regulation (EU) 10/2011 for short-term food contact



Item number	Material	Colour	Material thickness	Length x width	Design	Multipack	Pieces per carton	Cartons per pallet
703045-b	Polyethylene (PE)	Blue	24 µm	120 x 75 cm	Blocked	50 pieces	500 pieces	48
703045-w	Polyethylene (PE)	White	24 µm	120 x 75 cm	Blocked	50 pieces	500 pieces	40



medi-inn® CLASSIC PE Apron, PLUS

Non-sterile PE apron, for single use only

- Made of polyethylene (PE)
- Hammered
- Liquid-proof
- Neck cut-out and waistbands
- Latex-free
- Non-sterile
- Category I personal protective equipment in accordance with Regulation (EU) 2016/425



Item number	Surface	Colour	Material thickness	Length x width	Design	Multipack	Pieces per carton	Cartons per pallet
N12989	Hammered	White	16 µm	125 x 80 cm	Individually folded in dispenser box	100 pieces	1000 pieces	18
N13461	Hammered	White	16 µm	140 x 80 cm	Individually folded in dispenser box	80 pieces	800 pieces	18

medi-inn® CLASSIC PE Apron, PLUS

Non-sterile PE apron, for single use only, blocked



- Made of polyethylene (PE)
- Hammered
- Liquid-proof
- Neck cut-out and waistbands
- Latex-free
- Non-sterile
- Category I personal protective equipment in accordance with Regulation (EU) 2016/425



Item number	Surface	Colour	Material thickness	Length x width	Design	Multipack	Pieces per carton	Cartons per pallet
N11318-w	Hammered	White	17 µm	140 x 75 cm	50 pieces blocked, with tear-off perforation	50 pieces	500 pieces	40

medi-inn® CLASSIC Stainless Steel Holder for Aprons in Box



Wall bracket for aprons in boxes

- Stainless steel
- For wall mounting
- Suitable for storage and dispensing of aprons in the box
- Suitable for BODY PRODUCTS medi-inn® CLASSIC PE apron, PLUS, item N12989 & N13461



Item number	Material	Colour	Dimensions	Single pack	Pieces per carton	Cartons per pallet
N13504-STD	Stainless steel	Silver matt	Height: 110 mm Width: 275 mm Depth: 125 mm	1 piece	6 pieces	27

medi-Inn® CLASSIC Stainless Steel Holder for Blocked Aprons



Adjustable wall bracket for blocked aprons

- Stainless steel
- For wall mounting
- Suitable for storage and dispensing of blocked aprons
- Universally suitable for all apron manufacturers




Item number	Material	Colour	Dimensions	Single pack	Pieces per carton	Cartons per pallet
N13505-STD	Stainless steel	Silver matt	Height: 65 mm Width: 360 mm Depth: 10 mm	1 piece	6 pieces	42

medi-inn® CLASSIC Face Mask Type II

Medical face masks made of soft and skin-friendly non-woven fabric, 3-layer with integrated nose clip



- Type II in accordance with EN 14683:2019+AC:2019
- Bacterial filtration performance $\geq 98\%$
- Breathing resistance less than 40 Pa / cm²
- Available in 4 colours
- Elastic ear loops
- Adjustable nose clip
- Fibreglass-free
- Latex-free
- Non-sterile
- Class I medical device in accordance with Regulation (EU) 2017/745



Item number	Material	Colour	Layers	Length x width	Contents of a box	Pieces per carton	Cartons per pallet
703056-b-1000	Polypropylene non-woven	Blue	3-layers	95 x 175 mm	50 pieces	1000 pieces	24
703056-g-1000	Polypropylene non-woven	Green	3-layers	95 x 175 mm	50 pieces	1000 pieces	24
703056-w	Polypropylene non-woven	White	3-layers	95 x 175 mm	50 pieces	500 pieces	48
703056-s	Polypropylene non-woven	Black	3-layers	95 x 175 mm	50 pieces	500 pieces	40



medi-inn® CLASSIC Face Mask Type IIR

Medical face masks made of soft and skin-friendly non-woven fabric, 3-layer, with integrated nose clip

- With 4 straps for tying or ear loops
- Bacterial filtration performance $\geq 98\%$
- Breathing resistance less than 60 Pa / cm²
- Resistant to penetration of liquid ≥ 16 kPa
- Type IIR in accordance with EN 14683:2019+AC:2019
- Adjustable nose clip
- Fibreglass-free
- Latex-free
- Non-sterile
- Class I medical device in accordance with Regulation (EU) 2017/745



Item number	Material	Colour	Design	Layers	Length x width	Contents of a box	Pieces per carton	Cartons per pallet
N12359-g	Polypropylene non-woven	Green	Ear loops	3-layers	95 x 175 mm	50 pieces	500 pieces	40
N18240-g	Polypropylene non-woven	Green	For tying	3-layers	95 x 175 mm	50 pieces	500 pieces	32



medi-inn® GREEN Face Mask Type IIR

Biodegradable medical face mask made of polylactide (PLA) non-woven fabric, 3-layer, with integrated nose clip

- Ear loops, elasticated straps covered with cotton
- The face mask was tested in accordance with GB/T 19277.1-2011 and achieved a biodegradability of up to 90% in 130 days
- Resistant to penetration of liquid ≥ 16 kPa
- Type IIR in accordance with EN 14683:2019+AC:2019
- Flexible nose clip
- Fibreglass-free
- Latex-free
- Non-sterile
- Class I medical device in accordance with Regulation (EU) 2017/745



Item number	Material	Colour	Layers	Length x width	Contents of a box	Pieces per carton	Cartons per pallet
N640539	Poly lactide (PLA) nonwoven	White	3-layers	95 x 175 mm	50 pieces	500 pieces	35



medi-inn® CLASSIC Respirator FFP2 NR / Type IIR

FFP2 respirator with elasticated head straps, 4-layer

- Total inward leakage ≤ 8%
- Maximum flow of the filter medium ≤ 6%
- Type IIR: EN 14683:2019+AC:2019
- FFP2 NR: EN 149:2001+A1:2009
- Inhalation resistance 30 l/min ≤ 0.7 mbar
Inhalation resistance 95 l/min ≤ 2.4 mbar
Inhalation resistance 160 l/min ≤ 3.0 mbar
- Adjustable nose clip
- Head loops
- Latex-free
- Individually packed
- Non-sterile
- Class I medical device in accordance with Regulation (EU) 2017/745
- Category III personal protective equipment in accordance with Regulation (EU) 2016/425



Item number	Material	Colour	Layers	Length x width	Contents of a box	Pieces per carton	Cartons per pallet
N640574	Polypropylene non-woven	White	4-layers	160 x 107 mm	50 pieces	500 pieces	14



medi-inn® CLASSIC Respirator FFP2 NR

FFP2 respirator with elasticated ear loops, 4-layer

- Maximum flow of the filter medium $\leq 6\%$
- FFP2 NR: EN 149:2001+A1:2009
- Inhalation resistance 30 l/min ≤ 0.7 mbar
Inhalation resistance 95 l/min ≤ 2.4 mbar
Inhalation resistance 160 l/min ≤ 3.0 mbar
- Adjustable nose clip
- Elastic ear loops
- Latex-free
- Individually packed
- Non-sterile
- Category III personal protective equipment in accordance with Regulation (EU) 2016/425



Item number	Material	Colour	Layers	Length x width	Contents of a bag	Pieces per carton	Cartons per pallet
N18778	Polypropylene non-woven	White	4-layers	155 x 109 mm	5 pieces	500 pieces	18



medi-inn® CLASSIC Shoe Covers

Made of Cast Polyethylene with Elastic Band

CPE shoe covers

- Shoe covers made of cast polyethylene fabric with elastic band, hammered surface
- Suitable for shoe sizes 38 to 47
- Single-use
- Non-sterile
- Class I medical device in accordance with Regulation (EU) 2017/745




Item number	Material	Colour	Material thickness	Length x width	Design	Multipack	Pieces per carton	Cartons per pallet
703021	Cast polyethylene (CPE)	Blue	20 µm	41 x 15 cm	Packed in batches of 20 pieces	100 pieces	2000 pieces	18

medi-inn® CLASSIC Anti-Slip Socks



Anti-slip patient socks with anti-slip grip dots

- Double-sided rubberized
- Anti-slip grip dots
- Material sock: Polyester (PES)
Material sole: Polyvinyl chloride (PVC)
- Latex-free
- Non sterile
- Packed in pairs in a bag
- Available in sizes S-XXXL



Item number	Material	Colour	Size	Weight	Multipack	Cartons per pallet
N13958-S	Polyester	● Red	S (shoe size 32-35)	22 g	48 Individual packs x 1 Pair	77 cartons (3,696 pairs)
N13958-M	Polyester	● Green	M (shoe size 36-38)	24 g	48 Individual packs x 1 Pair	72 cartons (3,456 pairs)
N13958-L	Polyester	● Blue	L (shoe size 39-42)	30 g	48 Individual packs x 1 Pair	56 cartons (2,688 pairs)
N13958-XL	Polyester	● Brown	XL (shoe size 43-44)	37 g	48 Individual packs x 1 Pair	56 cartons (2,688 pairs)
N13958-g-XL	Polyester	● Green	XL (shoe size 43-44)	37 g	48 Individual packs x 1 Pair	48 cartons (2,304 pairs)
N13958-XXL	Polyester	● Grey	XXL (shoe size 45-46)	43 g	48 Individual packs x 1 Pair	40 cartons (1,920 pairs)
N13958-XXXL	Polyester	● Black	XXXL (shoe size 47-49)	48 g	48 Individual packs x 1 Pair	40 cartons (1,920 pairs)

Sizes S and M are children's sizes



Ward Supplies

- Denture Containers
- Patient Belongings Bags
- Medication Administration
- Medication Dispensers
- Medication Cutters
- Urinal Bottles and Accessories
- Bedpans
- Kidney Dishes
- Emesis Bags
- Cosmetics



medi-inn® CLASSIC Dental Containers

Container with insert and screw cap for transporting

- Denture container with screw lid
- 3-part: Lid, sieve insert and container
- Autoclavable up to 121 °C
- Material: Polypropylene (PP)
- Non-sterile



Item number	Material	Colour	Dimensions (diameter / height)	Single pack	Pieces per carton	Cartons per pallet
200100	Polypropylene (PP)	White	Cup: 9.6 cm / 5.4 cm Insert (sieve): 9.2 cm / 5.4 cm Lid: 9.6 cm / 1.7 cm	1 piece	20 pieces	54

medi-inn® CLASSIC Patient Belongings Bag

For the safe and organized storage of personal patient belongings



- Practical carry handle for easy transport
- Large labeling area with pre-printed lines for name, department, address, and date
- Durable yet flexible material
- Made of low-density polyethylene (LDPE)




Item number	Material	Colour	Dimensions	Single pack	Pieces per carton	Cartons per pallet
N11361	Polyethylen (LDPE)	White with blue print	58 cm x 50 cm x 0.4 cm	250 pieces	250 pieces	52



medi-inn® CLASSIC Medication Dispenser STANDARD

Medication dispenser with 4 compartments and sliding blister cover

- With transparent sliding blister
- Compartments for morning, noon, evening, night
- Printed with: "Morgen, Mittag, Abend und Nacht" / "morning - noon - evening - night"
- Embossed '1 2 3 4' incl. Braille
- Material: Polypropylene (PP), sliding blister made of polyvinyl chloride (PVC)
- Free from phthalates
- Non-sterile



Item number	Materials	Colour	Dimensions (length x width x height)	Multipack	Pieces per carton	Cartons per pallet
N17634-b-STD	Polypropylene (PP), sliding blister made of polyvinyl chloride (PVC)	Blue	Dispenser 22 cm x 6 cm x 1.7 cm Tablet compartment 4 cm x 2.5 cm x 1.6 cm	40 pieces	200 pieces	18
N17634-w-STD	Polypropylene (PP), sliding blister made of polyvinyl chloride (PVC)	White	Dispenser 22 cm x 6 cm x 1.7 cm Tablet compartment 4 cm x 2.5 cm x 1.6 cm	40 pieces	200 pieces	18



medi-inn® CLASSIC Medication Dispenser XL SPECIAL

Medication dispenser with 4 compartments and sliding blister cover

- With transparent sliding blister
- Compartments for morning, noon, evening, night
- Printed with: "Morgen, Mittag, Abend und Nacht" / "morning - noon - evening - night"
- Embossed '1 2 3 4' incl. Braille
- Material: Polypropylene (PP), sliding blister made of polyvinyl chloride (PVC)
- Free from phthalates
- Non-sterile



Item number	Materials	Colour	Dimensions (length x width x height)	Single pack	Pieces per carton	Cartons per pallet
N17672-b	Polypropylene (PP), sliding blister made of polyvinyl chloride (PVC)	Blue	Dispenser 22.6 cm x 8 cm x 1.5 cm Tablet compartment 5.8 cm x 4.7 cm x 1.7 cm	1 piece	190 pieces	18
N17672-w	Polypropylene (PP), sliding blister made of polyvinyl chloride (PVC)	White	Dispenser 22.6 cm x 8 cm x 1.5 cm Tablet compartment 5.8 cm x 4.7 cm x 1.7 cm	1 piece	190 pieces	18

medi-inn® CLASSIC M Cup ECO

Cups for collecting and distributing tablets or liquid medication



- Graduated from 5 ml to 25 ml in 1 ml increments
- Volume: 30 ml
- Intended for single use
- Material: Polypropylene (PP)
- Non-sterile
- Class Im medical device in accordance with Directive 93/42 EEC



CE0197



Item number	Material	Colour	Dimensions	Multipack	Pieces per carton	Cartons per pallet
N19032-b	Polypropylene (PP)	Blue	Diameter top: 38 mm Diameter bottom: 27 mm Height: 46 mm	80 pieces	4800 pieces	26
N19032-g	Polypropylene (PP)	Green	Diameter top: 38 mm Diameter bottom: 27 mm Height: 46 mm	80 pieces	4800 pieces	26
N19032-n	Polypropylene (PP)	Transparent	Diameter top: 38 mm Diameter bottom: 27 mm Height: 46 mm	80 pieces	4800 pieces	26
N19032-r	Polypropylene (PP)	Red	Diameter top: 38 mm Diameter bottom: 27 mm Height: 46 mm	80 pieces	4800 pieces	26
N19032-y	Polypropylene (PP)	Yellow	Diameter top: 38 mm Diameter bottom: 27 mm Height: 46 mm	80 pieces	4800 pieces	26

medi-inn® CLASSIC Lid M Cup ECO

Lid for closing



- Material: Polypropylene (PP)
- Non-sterile



Item number	Material	Colour	Dimensions	Multipack	Pieces per carton	Cartons per pallet
N19033	Polypropylene (PP)	White	Diameter top: 40 mm Height: 4 mm	500 pieces	4000 pieces	26

medi-inn® CLASSIC Medication Tray



For storing and transporting medication dispensers and/or medicine cups

- Holds up to 19 medication dispensers
- Holds up to 19 medicine cups
- Non-sterile



Item number	Material	Colour	Dimensions (length x width x height)	Single pack	Pieces per carton	Cartons per pallet
100700	Polypropylene (PP)	White	44 cm x 34 cm x 1.5 cm	1 piece	30 pieces	8

medi-inn® CLASSIC Tablet Crusher with Integrated Tablet Compartment



For quick and easy crushing of tablets

- Tablet crusher with storage compartment for tablets
- Collecting compartment for crushed tablets or tablet powder
- Material: Polypropylene (PP)
- Weight: 52 g
- Dishwasher safe up to 90°C
- Reusable
- Individually packed in PE bag and box
- Non-sterile



Item number	Materials	Colour	Dimensions	Single pack	Pieces per carton	Cartons per pallet
100122-b	Polypropylene (PP), polystyrene (PS)	Blue	Width: 4.5 cm Length: 5.2 cm	1 piece	100 pieces	32



medi-inn® CLASSIC Weekly Dispenser 7 Days

Weekly dispenser for storing medication

- Days of the week printed on the compartments
- Daily compartments printed with: Morning, Noon, Evening, Night
- Material: Acrylonitrile butadiene styrene copolymers (ABS), polystyrene (PS), polypropylene (PP)
- Free from phthalates
- Non-sterile
- Individually packed in folding boxes



Item number	Materials	Colour	Dimensions (length x width x height)	Single pack	Pieces per carton	Cartons per pallet
N11484-w	Acrylonitrile butadiene styrene copolymers (ABS), polystyrene (PS), polypropylene (PP)	White	Dispenser 11.4 cm x 4.3 cm x 12.8 cm Tablet compartment: 2.5 cm x 3.0 cm x 1.5 cm	1 piece	40 pieces	18

medi-inn® CLASSIC PP Urinal Bottle for Men



Male urinal bottle made of polypropylene, translucent

- Bottle autoclavable up to 134°C (without lid)
- Translucent
- Square design
- With blue lid
- Capacity 1,000 ml
- Material urinal bottle: polypropylene (PP)
Material lid: polypropylene (PP)
- Non-sterile
- Class I medical device in accordance with Regulation (EU) 2017/745

CE
MD

Item number	Material	Colour	Dimensions	Single pack	Pieces per carton	Cartons per pallet
100481	Polypropylene (PP)	Bottle: Translucent Lid: Blue	Length: 31.5 cm Base: 9 cm x 9.5 cm x 9.5 cm Bottleneck: Ø 5.5 cm Lid: Ø 5.9 cm	1 piece	25 pieces	10

medi-inn® CLASSIC PC Urinal Bottle for Men



Male urinal bottle made of polycarbonate, crystal clear

- Bottle autoclavable up to 134°C (without lid)
- Crystal clear / transparent
- Square design
- With blue lid
- Capacity 1,000 ml
- Material urine bottle: polycarbonate (PC)
Material lid: polypropylene (PP)
- Non-sterile
- Class I medical device in accordance with Regulation (EU) 2017/745

CE
MD

Item number	Materials	Colour	Dimensions	Single pack	Pieces per carton	Cartons per pallet
100483	Polycarbonate (PC) Polypropylene (PP)	Bottle: Crystal clear Lid: Blue	Length: 31.5 cm Base: 9 cm x 9.5 cm x 9.5 cm Bottleneck: Ø 5.5 cm Lid: Ø 5.9 cm	1 piece	25 pieces	10

medi-inn® CLASSIC PP Urinal Bottle for Women



Female urinal bottle made of polypropylene, translucent

- Bottle autoclavable up to 134°C (without lid)
- Translucent
- Square design
- With blue lid
- Capacity 1,000 ml
- Material urine bottle: polypropylene (PP)
Material lid: polypropylene (PP)
- Non-sterile
- Class I medical device in accordance with Regulation (EU) 2017/745



Item number	Material	Colour	Dimensions	Single pack	Pieces per carton	Cartons per pallet
N11516	Polypropylene (PP)	Bottle: Translucent Lid: White	Length: 26.5 cm Base: 14.8 cm x 6.8 cm Bottleneck (oval): 8.3 cm x 4.7 cm Bottleneck (oval): 9.5 cm x 6.0 cm	1 piece	48 pieces	4

medi-inn® CLASSIC Urinal Bottle Holder



Urinal bottle holder with hook for bed attachment and lid secured with a chain

- Only suitable for male urinal bottles
- Lid with chain
- Holder with hook for bed frame attachment
- Plastic-coated 4 mm wire
- For bed frames with a thickness of up to 5.5 cm
- Non-sterile



Item number	Colour	Dimensions	Single pack	Pieces per carton	Cartons per pallet
100484	White	Length: 32.5 cm Width: 15.5 cm	1 piece	25 pieces	18

medi-inn® CLASSIC Replacement Lid for Male Urinal Bottles



Lid for sealing urine bottles

- Spare part for urine bottles
- Material: polypropylene (PP)
- Non-sterile



Item number	Material	Colour	Dimensions	Multipack	Pieces per carton	Cartons per pallet
100489-b	Polypropylene (PP)	Blue	Ø 5.9 cm	25 pieces	250 pieces	55

medi-inn® CLASSIC Cleaning Brush



Brush for cleaning

- Flexible wire brush with plastic bristles
- Long handle for thorough cleaning
- Non-sterile



Item number	Colour	Dimensions	Multipack	Pieces per carton	Cartons per pallet
100490	Brush: White Handle: Blue	Total length: 37 cm Brush length: 19 cm Diameter top: 5.5 cm	25 pieces	250 pieces	4



medi-inn® CLASSIC Bedpan

Polypropylene (PP) bedpan with handle and lid, reusable

- Container for bedridden people to for urination and/or defecation
- Reusable
- Material: polypropylene (PP)
- With handle and lid
- Maximum load: 150 kg
- Non-sterile
- Class I medical device in accordance with Regulation (EU) 2017/745



Item number	Material	Colour	Dimensions	Multipack	Pieces per carton	Cartons per pallet
N11470-y	Polypropylene (PP)	Yellow	Length with handle: 430 mm Inner diameter: 190 mm Outer diameter: 310 mm Height: 100 mm	16 pieces	16 pieces	8

medi-inn® CLASSIC Disposable Kidney Dish



Disposable cardboard kidney dish

- For body fluids or small medical waste in healthcare settings
- Material: 100% recycled paper fibre
- Resistant to liquids for 24 hours:
Distilled water
Isotonic saline solution
Hydrogen peroxide 3%
Urea solution 3%
- Inner packed in polybags of 50 pieces
- Non-sterile
- Class I medical device in accordance with Regulation (EU) 2017/745



Item number	Material	Colour	Dimensions	Multipack	Pieces per carton	Cartons per pallet
N15000-n	Paper fibre	Brown-grey	Length: 255 ± 2 mm Width: 150 ± 2 mm Height: 45 ± 2 mm Capacity: 700 ml	50 pieces	300 pieces	18

medi-inn® CLASSIC Stainless Steel Kidney Dish



Stainless steel kidney dish, reusable

- Versatile application options
- Easy to clean
- Material: stainless steel
- RDG suitable and autoclavable up to 134°C
- Individually packed
- Non-sterile
- Class I medical device in accordance with Regulation (EU) 2017/745




Item number	Material	Colour	Dimensions	Multipack	Pieces per carton	Cartons per pallet
N19009	Stainless steel	Silver	Length: 265 ± 5 mm Width: 150 ± 3 mm Height: 38 ± 1 mm Capacity: 800 ml	5 pieces	50 pieces	18



medi-inn® CLASSIC Collection Bag Plus

Collection bags protect the user and surroundings from vomit and body fluids

- Tear-resistant polyethylene (PE)
- Protects the environment and users against contamination from vomit
- Easy to close
- Free from latex, DEHP and PVC
- Inner packed in polybags of 25 pieces
- Non-sterile



Item number	Material	Colour	Dimensions	Multipack	Pieces per carton	Cartons per pallet
N10345	Polyethylene (PE)	White	Closing ring: Ø 12.5 cm Length with bag: 36.5 cm Volume: approx. 4 litres	25 pieces	400 pieces	18

medi-inn® CLASSIC Comb

Comb for personal hygiene



- Lightweight and easy to carry
- Antistatic
- Made of plastic (polypropylene PP)



Item number	Material	Colour	Dimensions	Multipack	Pieces per carton	Cartons per pallet
N11034	Polypropylene (PP)	Black	Length side wide teeth: 5.8 cm Length side fine teeth: 5.5 cm Length (total): 12.5 cm	50 pieces	1000 pieces	42

medi-inn® CLASSIC Comb

Comb for personal hygiene



- Lightweight and easy to carry
- Antistatic
- Made of plastic (polypropylene PP)



Item number	Material	Colour	Dimensions	Multipack	Pieces per carton	Cartons per pallet
N11035	Polypropylene (PP)	Black	Length side wide teeth: 8.3 cm Length side fine teeth: 8.0 cm Length (total): 17.5 cm	50 pieces	1000 pieces	36



medi-inn® CLASSIC Cosmetic Kit

10-piece personal care set in a resealable, transparent pouch with zipper closure

- Resealable, transparent bag
- Contents:
 - Fine soap - 10 g
 - Shower gel - 30 ml
 - Hair shampoo - 30 ml
 - Body lotion - 30 ml
 - Terry towel - approx. 30 cm x 30 cm
 - Shower cap with elastic band
 - Toothbrush - individually packed in clear bag
 - Toothpaste - 6 g
 - Comb
 - Nail file



Item number	Colour	Pouch dimensions (mm)	Single pack	Pieces per carton	Cartons per pallet
N17686-BPR	Transparent	Length: 18.5 cm Width: 5.5 cm Height: 13.0 cm	1 case	20 cases	28

The cosmetic formulations are labelled and registered in accordance with Regulation (EC) No. 1223/2009 on cosmetic products and comply with all applicable requirements. The cosmetic products included are nanoparticle-free.



04 Diagnostics Aids

- Stethoscopes
- Blood Pressure Monitor
- Thermometers
- Urine Specimen Container
- Tourniquets
- Telemetry Pouches
- Wooden Tongue Depressors
- Diagnostic Examination Lights



medi-inn® CLASSIC Rappaport Stethoscope with Dual Head

The stethoscope can be adjusted in five different ways and offers maximum flexibility for various diagnostic purposes.

- Dual-tube system and swivelling double head
- Spring-loaded ear tubes for a comfortable fit and secure fit
- Set with 3 additional chest pieces in various sizes
- Includes 2 replacement diaphragms (small and large) and 2 pairs of replacement ear tips
- Can be adjusted in five different ways for different diagnostic requirements
- Non-sterile
- Class I medical device in accordance with Regulation (EU) 2017/745



Also included in the delivery:
3 breast pieces,
2 replacement membranes,
2 pairs of replacement ear buds



Item number	Length	Tube length	Ear tube length	Diaphragm diameter	Colour	Multipack	Pieces per carton
710039-s-5	735 mm	540 mm	160 mm	45.5 / 31 mm	Black	5 pieces	50 pieces

medi-inn® CLASSIC Stethoscope Single Head



Robust flat-head stethoscope for differentiated examinations

- Durable, noise-insulating PVC tubing
- Y-shaped tubing
- Flat chest piece for precise diagnostics
- Materials: aluminium and polyvinyl chloride (PVC)
- Non-sterile
- Class I medical device in accordance with Regulation (EU) 2017/745



Item number	Length	Tube length	Ear tube length	Diaphragm diameter	Colour	Multipack	Pieces per carton
710040-b-10	750 mm	600 mm	160 mm	46 mm	Blue	10 pieces	100 pieces
710040-g-10	750 mm	600 mm	160 mm	46 mm	Green	10 pieces	100 pieces
710040-r-10	750 mm	600 mm	160 mm	46 mm	Red	10 pieces	100 pieces
710040-s-10	750 mm	600 mm	160 mm	46 mm	Black	10 pieces	100 pieces
710040-y-10	750 mm	600 mm	160 mm	46 mm	Yellow	10 pieces	100 pieces

medi-inn® CLASSIC Duale Head Stethoscope



Versatile dual head stethoscope with optimum fit

- Nickel-free
- Black Y-shaped tubing
- Rotating dual head allows selection of desired stethoscope side
- Materials: aluminium and polyvinyl chloride (PVC)
- Non-sterile
- Class I medical device in accordance with Regulation (EU) 2017/745



Item number	Length	Tube length	Ear tube length	Diaphragm diameter	Colour	Multipack	Pieces per carton
710045-s-10	780 mm	615 mm	160 mm	46.5 mm 29.8 mm	Black	10 pieces	100 pieces



medi-inn® CLASSIC Sphygmomanometer Aneroid Premium, 1- or 2-Tube Version

Blood pressure measurement in healthcare facilities, at home, at work and in educational institutions

- Shock- and overpressure-resistant manometer
- Pressure gauge in metal housing
- Black nylon cuff, suitable for arm circumferences from 25 to 40 cm
- Supplied with practical zip pouch
- Non-sterile
- Class I medical device with measuring function in accordance with directive 93/42/EEC



CE0197



Item number	System	Measuring range	Measuring resolution	Scale diameter	Cuff size	Single pack	Pieces per carton
N13542	1-hose	0 - 300 mmHg	2 mmHg	60 mm	25 - 40 cm	1 piece	20 pieces
N11564	2-hose	0 - 300 mmHg	2 mmHg	60 mm	25 - 40 cm	1 piece	20 pieces



medi-inn® CLASSIC Blood Pressure Monitor for Upper Arm or Wrist

Fully automatic digital blood pressure monitors, including storage pouch and suitable AAA batteries

- Oscillometric measurement
- Automatic inflation
- Auto shut-off function
- Arrhythmia detection
- Memory for 2 users with 150 memory slots
- Powered by 2x AAA (wrist model) or 3x AAA batteries (upper arm model), batteries included; optional 5 V power adapter (not included)
- Non-sterile
- Class IIa medical device in accordance with Regulation (EU) 2017/745



CE0123 MD


Item number	System	Dimensions	Screen	Measuring range	Measuring accuracy	Cuff size	Single pack	Pieces per carton
N640502	Upper arm	136 x 94.5 x 57 mm Hose length: 600 mm	50 x 62 mm (LCD)	Systolic 60 - 260 mmHg Diastolic 30 - 200 mmHg Pulse 30 - 180 min ⁻¹	Pressure ± 3 mm Pulse ± 5%	22 - 26 cm	1 piece	24 pieces
N640504	Wrist	84 x 62 x 23.5 mm	32 x 21.5 mm (LCD)	Systolic 60 - 260 mmHg Diastolic 30 - 200 mmHg Pulse 30 - 180 min ⁻¹	Pressure ± 3 mm Pulse ± 5%	13.5 - 21.5 cm	1 piece	48 pieces



medi-inn® CLASSIC Clinical Thermometer – Rigid Tip

- Measuring range 32 - 42.9°C
- Measuring tip made of stainless steel
- Quick measurement: signal after approx. 10 seconds
- Last value memory
- IP protection class 22
- Individually packed
- Non-sterile
- Class IIa medical device in accordance with Regulation (EU) 2017/745



CE0123 MD 

Item number	System	Size	Tip diameter	Measuring range	Materials	Multipack	Pieces per carton
N11038-10	rigid	12.3 cm x 0.9 cm x 1.8 cm	3 mm	32 - 42.9°C	ABS, stainless steel	10 pieces	100 pieces

medi-inn® CLASSIC Clinical Thermometer – Flexible Tip



- Measurement range 32 - 42.9°C
- Measuring tip made of stainless steel
- Quick measurement: signal after approx. 10 seconds
- Last value memory
- IP protection class 27
- Individually packed
- Non-sterile
- Class IIa medical device in accordance with Regulation (EU) 2017/745



CE0123 MD 

Item number	System	Size	Tip diameter	Measuring range	Materials	Multipack	Pieces per carton
N14862-10	flexible	12.4 cm x 1.2 cm x 2.4 cm	3 mm	32 - 42.9°C	ABS, TPE (tip) stainless steel	10 pieces	100 pieces



medi-inn® CLASSIC Urine Specimen Container 125 ml, with Lid and Red Spout

Storage of urine samples for in-vitro diagnostics

- Sealable container for collecting body fluids
- Capacity: 125 ml
- Graduation in 25 ml increments
- Single use
- Non-sterile
- In vitro diagnostics medical device class A in accordance with Regulation (EU) 2017/746



Item number	Material	Colour	Dimensions	Multipack	Pieces per carton
100402-r	Polypropylene (PP)	Container: Transparent Lid: Red	Container Ø top: 63 mm Ø bottom: 48 mm Height: 72 mm Lid Height: 13 mm Ø: 65 mm	500 pieces	500 pieces



medi-inn® CLASSIC Urine Specimen Container 100 ml, with Green Screw Cap

- Sealable container for collecting body fluids
- Capacity: 100 ml
- Graduation in 10 ml increments
- Intended for single use
- Non-sterile
- In vitro diagnostics medical device class A in accordance with Regulation (EU) 2017/746



Item number	Material	Colour	Dimensions	Multipack	Pieces per carton
100451-g	Polypropylene (PP)	Container: Transparent Lid: Green	Container Ø top: 63 mm Ø bottom: 48 mm Height: 72 mm Lid Height: 13 mm Ø: 65 mm	5 x 100 pieces	500 pieces

medi-inn® CLASSIC Tourniquet



- Easy single-handed operation
- Individually packaged in folding box
- Non-sterile
- Class I medical device in accordance with Regulation (EU) 2017/745



Item number	Size	Closure	Colour	Band Materials	Closure Materials	Multipack	Pieces per carton
704003-b-10	45 x 2.5 cm	5.6 x 3.2 cm	Blue	Polyester, spandex	Acrylonitrile butadiene styrene (ABS)	10 pieces	100 pieces

medi-inn® CLASSIC Telemetry Pouch



For the secure storage and carrying of monitoring devices for long-term ECG, long-term blood pressure and similar monitoring devices

- Velcro closure at the top
- 4 straps for fixation on the body
- Single use
- Non-sterile



Item number	Material	Material thickness	Colour	Dimensions	Multipack	Pieces per carton
N15655-w	polypropylene (PP) Spunbond	13 my	White - Transparent	Outer dimensions of pouch: 21 x 16 cm (L x W) Inner dimensions of pouch: 17 x 14 cm (L x W) Straps: 102 x 1 cm (L x W)	50 pieces	300 pieces

medi-inn® CLASSIC Wooden Tongue Depressors



Wooden mouth spatulas for diagnostic use

- Individually removable diagnostic spatulas made of birch wood with rounded edges
- For inspection of the throat
- Single use
- Non-sterile (not individually wrapped)
- Class I medical device in accordance with Regulation (EU) 2017/745



Item number	Material	Colour	Dimensions	Multipack	Pieces per carton
NT1924	Birch wood	Brown	Length: 150 mm Width: 18 mm Height: 1.6 mm	100 pieces	300 pieces

medi-inn® CLASSIC Diagnostic Penlight



- For diagnostics of the oral cavity, throat and skin surface
- Not suitable for direct eye illumination
- Clip for easy attachment
- Individually packed in blisters
- 2 AAA batteries included in the delivery
- Non-sterile



Item number	Length	Brightness	Colour	Material	Multipack	Pieces per carton
704001	12.5 cm	8 lumen	White	PET	50 pieces	500 pieces



05

Nursing Care

- Waterless Body Washing
- Wash Gloves
- Hair Removal
- Oral Care
- Bibs
- Feeding Cups

medi-inn® CLASSIC Moist Wash Gloves



Moist wash gloves in resealable polybag

- Soft and skin-friendly
- pH-neutral
- Alcohol-free, EDTA-free, paraben-free, nanoparticle-free
- Non-sterile
- Latex-free
- According to Cosmetics Regulation (EC) No. 1223/2009



Item number	Colour	Grammage / weight	Dimensions	Multipack	Pieces per carton	Cartons per pallet
701142-ALOE	White	80 g / m ²	Length: 23.5 cm Width: 14.5 cm	8 pieces	192 pieces	45

medi-inn® CLASSIC Cleansing Foam, 500 ml



The wash- and cleansing foam gently cleanses the skin without the use of soap additives or water

- Soft and skin-friendly
- Alternative to soap and water
- Suitable for daily use
- Effective cleansing components
- pH-neutral
- Nanoparticle-free
- Contents: 500 ml
- Non-sterile
- According to Cosmetics Regulation (EC) No. 1223/2009



Item number	Single pack	Pieces per carton	Cartons per pallet
N10957	1 piece	12 pieces	56



medi-inn® CLASSIC Moist Care Wipes (Disposable)

Disposable moist care wipes in dispenser packs for skin cleansing. Individually removable, resealable

- Soft and skin-friendly
- Dermatologically tested
- Practical for use at home and on the go
- Resealable soft pack with white plastic lid
- Fragrance- and paraben-free, pH-neutral
- Nanoparticle-free
- Non-sterile
- According to Cosmetics Regulation (EC) No. 1223/2009
- 80 wipes per pack



Item number	Material	Colour	Grammage / weight	Dimensions	Contents per soft pack	Pieces per carton	Cartons per pallet
N13457	Non-woven spunlace	White	45 g / m ²	Width: 18 ± 1 mm Length: 20 ± 1 mm	80 pieces	1280 pieces	28
N14235	Non-woven fabric	White	40 g / m ²	Width: 20 ± 1 mm Length: 27 ± 1 mm	80 pieces	960 pieces	30

medi-inn® CLASSIC Molton Wash Glove



Soft and tear-resistant wash glove made of polyester for single use

- Practical vacuum packaging
- Textile-like white molton, ultra-soft
- Stable edge seam, tear-resistant when wet
- Ideal for use with care products, does not fluff
- The pocket shape ensures comfortable, ergonomic handling
- Latex-free
- Non-sterile



Item number	Material	Colour	Grammage / weight	Dimensions	Multipack	Pieces per carton	Cartons per pallet
701158	Polyester non-woven fabric	White	75 g / m ²	Length: 22 ± 0.5 mm Width: 15 ± 0.5 mm	50 pieces	1000 pieces	18

medi-inn® CLASSIC Airlaid PE Wash Glove



Soft and tear-resistant wash glove made of airlaid with polyethylene (PE) for single use
Internal coating for single use

- Soft, absorbent skin-friendly airlaid top layer
- Low-linting, tear-resistant, highly absorbent
- The PE inner coating protects the hands during washing
- Non-sterile
- Packed in PE bag



Item number	Materials	Colour	Grammage / weight	Dimensions	Multipack	Pieces per carton	Cartons per pallet
701161	Cellulose (airlaid) + PE inner coating	Light blue	Airlaid: 50 g / m ² PE: 17 g / m ²	Length: 22 ± 0.5 mm Width: 16 ± 0.5 mm	50 pieces	1000 pieces	24



medi-inn® CLASSIC Wash Glove Economy

Soft polyester wash glove with high seam strength for single use only

- Gentle on the skin and lint-free
- Hygienic, highly absorbent and tear-resistant
- Soft, fluffy nonwoven fabric
- Non-sterile
- Practical vacuum packaging



Item number	Material	Colour	Grammage / weight	Dimensions	Multipack	Pieces per carton	Cartons per pallet
701172	Polyester non-woven fabric	White	60 g / m ²	Length: 22 ± 0.5 mm Width: 15 ± 0.6 cm	50 pieces	2000 pieces	18

medi-inn® CLASSIC Disposable Razor with Blade Guard, Single-Blade or Double-Blade



Single-blade or double-blade disposable razor for hygienic hair removal

- Stainless
- Single-use only
- Non-sterile
- Packed in a practical dispenser box with 100 pieces for individual removal



Item number	Design	Materials	Colour	Dimensions	Multipack	Pieces per carton	Cartons per pallet
704027	Double-blade	Handle: Polystyrene Blade: Stainless steel Protective cover: Polypropylene	Blue	Length: 10 cm Width (head): 4 cm	100 pieces	1000 pieces	15
704028	Single-blade	Handle: Polystyrene Blade: Stainless steel Protective cover: Polypropylene	Blue	Length: 10 cm Width (head): 4 cm	100 pieces	1000 pieces	15



medi-inn® CLASSIC Disposable Toothbrush with Toothpaste

Disposable toothbrush with transparent toothpaste

- Toothpaste is activated by contact with water
- Rounded edges and flexible bristles
- Blue handle
- Single-use only
- Non-sterile
- Individually hygienically packed
- 100 pieces in a practical dispenser box



Item number	Material	Colour	Weight	Dimensions	Multipack	Pieces per carton	Cartons per pallet
706200	Polypropylene (PP)	Blue with white bristles	5.0 g / piece	Length: 152 mm Width: 11 mm Brush: 13 mm long	100 pieces	1000 pieces	20

medi-inn® CLASSIC Oral Care Swabs



Sugar-free oral care swabs with foam head

- Sugar-free
- Soft foam for gentle oral hygiene and thorough particle removal
- For single use only
- Non-sterile
- Individually hygienically packed



Item number	Materials	Colour	Flavour	Dimensions	Multipack	Pieces per carton	Cartons per pallet
N12378-l	Swabs: Paper Head: Polyurethane foam	Foam head: Pink	Lemon	Length: 106 mm ± 4 mm Foam head: Ø: 18 mm ± 1 mm Length: 20 mm ± 0.5 mm	250 pieces	1000 pieces	21
N12378-m	Swabs: Paper Head: Polyurethane foam	Foam head: Pink	Mint	Length: 106 mm ± 4 mm Foam head: Ø: 18 mm ± 1 mm Length: 20 mm ± 0.5 mm	250 pieces	1000 pieces	21
N12378-n	Swabs: Paper Head: Polyurethane foam	Foam head: Green	Neutral	Length: 106 mm ± 4 mm Foam head: Ø: 18 mm ± 1 mm Length: 20 mm ± 0.5 mm	250 pieces	1000 pieces	21

medi-inn® CLASSIC Tooth Cup

Mug for use when brushing teeth



- White cup for storing toothbrushes and toothpaste
- For rinsing the mouth after brushing the teeth
- Dishwasher and microwave safe
- 250 ml graduated
- Non-sterile



Item number	Material	Colour	Dimensions	Multipack	Pieces per carton	Cartons per pallet
N17786	Polypropylene (PP)	White	Height: 94 mm Ø top: 72 mm Ø bottom: 58 mm	10 pieces	250 pieces	18

medi-inn® CLASSIC Bib

Disposable bib with crumb pocket, liquid-repellent



- Practical crumb pocket at the front
- Ideal for adults, seniors and children
- Protects clothing from soiling
- White tissue paper with blue PE backing
- Single-use
- Non-sterile



Item number	Material	Colour	Weight	Dimensions	Multipack	Pieces per carton	Cartons per pallet
730365-b	Tissue paper	Blue	Tissue: 18 g / m ² PE film: 11 g / m ²	Length: 64 cm Width: 37 cm	50 pieces	500 pieces	54
730365-g	Tissue paper	Green	Tissue: 18 g / m ² PE film: 11 g / m ²	Length: 64 cm Width: 37 cm	50 pieces	500 pieces	54



medi-inn® CLASSIC Feeding Cup – Base

Cup base for easier fluid intake

- Cup base with max. filling capacity: approx. 250 ml
- Dishwasher-safe up to 70 °C
- Compatible lids with different drinking openings available: 706302-t and 706304-t
- Non-sterile
- Class I medical device in accordance with Regulation (EU) 2017/745
- Suitable for food contact in accordance with Regulation (EC) 1935/2004



Item number	Material	Colour	Dimensions	Multipack	Pieces per carton	Cartons per pallet
706300-t	Polypropylene (PP)	Transparent	Height: 94 mm Ø: 70 mm	250 pieces	250 pieces	18



medi-inn® CLASSIC Feeding Cup Lid Ø 7 cm x 2 cm, Small Drinking Opening

For easier fluid intake when used with the feeding cup base 706300-t

- Drinking opening: 4 x 4 mm
- Lid for use with feeding cup base 706300-t
- Dishwasher-safe up to 70 °C
- Non-sterile
- Class I medical device in accordance with Regulation (EU) 2017/745
- Suitable for food contact in accordance with Regulation (EC) 1935/2004



Item number	Material	Colour	Dimensions	Multipack	Pieces per carton	Cartons per pallet
706302-t	Polypropylene (PP)	Transparent	Height: 20 mm Ø: 70 mm	250 pieces	250 pieces	18

medi-inn® CLASSIC Feeding Cup Lid Ø 7 cm x 2 cm, Large Drinking Opening

For easier fluid intake when used with the feeding cup base 706300-t



- Drinking opening: 12 x 10 mm
- Lid for use with feeding cup base 706300-t
- Dishwasher-safe up to 70 °C
- Non-sterile
- Class I medical device in accordance with Regulation (EU) 2017/745
- Suitable for food contact in accordance with Regulation (EC) 1935/2004



Item number	Material	Colour	Dimensions	Multipack	Pieces per carton	Cartons per pallet
706304-t	Polypropylene (PP)	Transparent	Height: 20 mm Ø: 70 mm	250 pieces	250 pieces	18




medi-inn® CLASSIC Feeding Cup Set, Transparent

Feeding cup set for easier fluid intake

- Set consisting of:
Lid with drinking opening Ø 4 mm
Lid with drinking opening Ø 12 mm
Cup base for approx. 250 ml
Packed in a polybag
- Practical patient set with two lids with differently sized drinking openings
- Dishwasher-safe up to 70 °C
- Non-sterile
- Class I medical device in accordance with Regulation (EU) 2017/745
- Suitable for food contact in accordance with Regulation (EC) 1935/2004



Item number	Material	Colour	Dimensions	Single pack	Sets per carton	Cartons per pallet
706350	Polypropylene (PP)	Transparent	Dimensions of individual components (706300-t, 706302-t and 706304-t) see previous page	1 set	250 pieces	18



06

Hygiene

- Examination Couch Rolls
- Patient Transfer Sheets
- Waste Disposal

medi-inn® CLASSIC Couch Protection Rolls

Disposable couch cover roll, 2-ply, white



- Versatile use: suitable not only for examination couches in medical practices, but also as a hygienic cover for massage tables, tattoo and cosmetic couches, as well as changing tables in childcare facilities
- Tear-resistant and highly absorbent
- Soft and skin-friendly
- 2-ply laminated (34 g/m²), perforated every 35 cm
- Total length: 47.25 metres per roll
- Available in widths: 39.5 cm, 49.5 cm, 55 cm, 59 cm
- Non-sterile
- Each roll is individually hygienically wrapped in film




Item number	Material	Material thickness	Colour	Dimensions	Single pack	Pieces per carton	Cartons per pallet
N10164	100% cellulose	17 g / m ²	White	Width: 39.5 cm Sheet length: 35 cm Total length: 47.25 m	1 piece	9 pieces	30
N10165	100% cellulose	17 g / m ²	White	Width: 49.5 cm Sheet length: 35 cm Total length: 47.25 m	1 piece	9 pieces	24
N10166	100% cellulose	17 g / m ²	White	Width: 55 cm Sheet length: 35 cm Total length: 47.25 m	1 piece	9 pieces	18
N12562	100% cellulose	17 g / m ²	White	Width: 59 cm Sheet length: 35 cm Total length: 47.25 m	1 piece	9 pieces	18



medi-inn® CLASSIC Patient Transfer Sheet

Thread-reinforced patient transfer sheet (8 threads) with liquid-tight underside, for single use

- Hygienic protection for stretchers and couches against moisture and contamination
- Tear-resistant and highly absorbent
- Soft and skin-friendly
- Liquid-proof PE backing on the underside
- 2-ply fabric, 24 g / m²
- Reinforced with 8 blue nylon threads
- Non-sterile
- Individually packed



Item number	Materials	Material thickness	Colour	Dimensions	Multipack	Pieces per carton	Cartons per pallet
730375	Paper + polyethylene (PE) film + nylon threads	Paper 24 g / m ² PE film 12 g / m ²	White	Length: 210 cm Width: 75 cm	25 pieces	100 pieces	42

medi-inn® CLASSIC Patient Transfer Sheet

Thread-reinforced patient transfer sheet (48 threads) with liquid-tight underside, for single use



- Hygienic protection for stretchers and couches against moisture and contamination
- Tear-resistant and highly absorbent
- Soft and skin-friendly
- Liquid-proof PE backing on the underside
- 2-ply fabric, 24 g / m²
- Reinforced with 48 blue nylon threads
- Non-sterile
- Individually packed



Item number	Materials	Material thickness	Colour	Dimensions	Multipack	Pieces per carton	Cartons per pallet
730376	Paper + polyethylene (PE) film + nylon threads	Paper 24 g / m ² PE film 12 g / m ²	White	Length: 210 cm Width: 75 cm	25 pieces	100 pieces	42

medi-inn® CLASSIC Patient Transfer Sheet

Thread-reinforced patient transfer sheet (8 threads) with liquid-tight underside, for single use

- Hygienic protection for stretchers and couches against moisture and contamination
- Tear-resistant and highly absorbent
- Soft and skin-friendly
- Liquid-proof PE backing on the underside
- 2-ply fabric, 24 g / m²
- Reinforced with 8 blue nylon threads
- Non-sterile
- Individually packed



Item number	Materials	Material thickness	Colour	Dimensions	Multipack	Pieces per carton	Cartons per pallet
741102	Paper + polyethylene (PE) film + nylon threads	Paper 24 g / m ² PE film 12 g / m ²	White	Length: 200 cm Width: 100 cm	25 pieces	100 pieces	28
741116	Paper + polyethylene (PE) film + nylon threads	Paper 24 g / m ² PE film 12 g / m ²	White	Length: 140 cm Width: 75 cm	25 pieces	100 pieces	60



medi-inn® CLASSIC Hygiene Bag Dispenser

Hygiene bag dispenser for plastic hygiene bags

- Easy wall mounting
- Including mounting material (2 screws, 2 plugs, 2 double-sided foam pads)
- Suitable for N11768 medi-inn CLASSIC hygiene bags and all common plastic hygiene bags (made of HDPE)
- Non-sterile
- Individually welded



Item number	Material	Colour	Dimensions	Single pack	Pieces per carton	Cartons per pallet
701124	Polypropylene (PP)	White	Length: 13.95 cm Height: 2.75 cm Width: 10 cm Dispenser opening: Ø 5.3 cm	1 piece	24 pieces	70
701128	Polypropylene (PP)	Silver (chrome-plated)	Length: 13.95 cm Height: 2.75 cm Width: 10 cm Dispenser opening: Ø 5.3 cm	1 piece	24 pieces	70

medi-inn® CLASSIC Hygiene Bag

For hygienic and convenient disposal of hygiene products



- Hygiene bags made of HDPE material (high-density polyethylene)
- Securely trap odours and liquids
- For the discreet disposal of sanitary towels, panty liners and tampons
- Non-sterile
- 30 pieces per pack



Item number	Material	Colour	Dimensions	Multipack	Pieces per carton	Cartons per pallet
N11768	Plastic (HDPE high-density polyethylene)	White	Width: 8 cm Height: 25 cm	30 pieces	1500 pieces	40



07

First Aid

- Emergency Blankets
- Instant Cold Packs
- Reusable Hot / Cold Packs



medi-inn® CLASSIC Emergency Blankets

The emergency blanket is designed to protect injured persons from cold or heat in emergency situations

- Protects against cold and heat
- Cold protection: silver side facing the body
- Heat protection: gold side facing the body
- Ultra-thin, lightweight, aluminium-coated film
- Non-sterile
- Individually packed in a polybag
- Class I medical device in accordance with Regulation (EU) 2017/745



Item number	Materials	Colour	Dimensions	Single pack	Pieces per carton	Cartons per pallet
851004	PET film, aluminium-coated	Gold / silver	Width: 160 cm Length: 210 cm Thickness: 12 µm	1 piece	200 pieces	27



medi-inn® CLASSIC Instant Cold Pack

Immediate help for bruises, sprains and haematomas, for external use only

- Ready to use in seconds
- Activation by squeezing the bag
- Provides even cold distribution
- Ideal for at home and on the go
- Non-sterile
- Individually packed
- Class IIa medical device in accordance with Regulation (EU) 2017/745



CE0197



Item number	Materials	Colour	Dimensions	Multipack	Pieces per carton	Cartons per pallet
N11970	Non-woven fabric with PE nylon	White-green	14 cm x 18 cm	25 pieces	25 pieces	48
N14843	PE nylon	White-green	12.7 cm x 15.3 cm	72 pieces	72 pieces	40
N14844	PE nylon	White-green	15.2 cm x 22.9 cm	36 pieces	36 pieces	40



medi-inn® CLASSIC Reusable Hot/Cold Pack

Reusable hot/cold pack for thermotherapy,
for external use, microwave-safe

- Can be used as a hot or cold pack
- Quick help for pain, strains and other complaints
- Freezer-safe and microwave safe
- Non-sterile
- Reusable
- Class I medical device in accordance with Regulation (EU) 2017/745



Not suitable for children under 3 years. Only to be used with adult supervision.

Item number	Materials	Colour	Dimensions	Weight	Multipack	Pieces per carton	Cartons per pallet
N18861	Outer film: Polyvinyl chloride (PVC) Gel: Water, glycerine, sodium carboxymethylcellulose	Blue	Length: 14 cm Width: 13 cm	140 g	10 pieces	60 pieces	47
N18869	Outer film: Polyvinyl chloride (PVC) Gel: Water, glycerine, sodium carboxymethylcellulose	Blue	Length: 29 cm Width: 12 cm	420 g	10 pieces	20 pieces	48
N18871	Outer film: Polyvinyl chloride (PVC) Gel: Water, glycerine, sodium carboxymethylcellulose	Blue	Length: 26 cm Width: 16 cm	550 g	10 pieces	20 pieces	48



Incontinence Care

- Pads
- Shaped Pads
- Incontinence Slips
- Pants
- Fixation Pants
- Protective Underpads
- Bed Protection

medi-inn® CLASSIC Pads

Light protection for active everyday life



- Shaped disposable pads
- Unisex fit
- Odour-binding and breathable
- Latex-free
- Colour coded packaging for easy absorbency level identification
- Non-sterile
- Class I medical device in accordance with Regulation (EU) 2017/745




Item number	Product	Colour coding for absorbency levels	Absorbency in accordance with ISO 11948-1	Multipack	Pieces per carton
N18552	medi-inn® CLASSIC Pads NORMAL	☹️	310 ml	28 pieces	168 pieces (6 x 28 pieces)
N18553	medi-inn® CLASSIC Pads PLUS	☹️☹️	410 ml	28 pieces	168 pieces (6 x 28 pieces)
N18555	medi-inn® CLASSIC Pads MAXI	☹️☹️☹️	820 ml	28 pieces	168 pieces (6 x 28 pieces)

medi-inn® CLASSIC Shaped Pads

Reliable protection with a secure fit



- Anatomically shaped disposable pads
- With wetness indicator
- Unisex fit
- Odour-binding and breathable
- Latex-free
- For secure fixation, we recommend combining them with our washable fixation pants
- Colour coded packaging for easy absorbency level identification
- Non-sterile
- Class I medical device in accordance with Regulation (EU) 2017/745




Item number	Product	Colour coding for absorbency levels	Absorbency in accordance with ISO 11948-1	Multipack	Pieces per carton
N18568	medi-inn® CLASSIC Shaped Pads NORMAL	👉👉👉	1,350 ml	30 pieces	120 pieces (4 x 30 pieces)
N18569	medi-inn® CLASSIC Shaped Pads PLUS	👉👉👉👉	1,700 ml	30 pieces	120 pieces (4 x 30 pieces)
N18570	medi-inn® CLASSIC Shaped Pads EXTRA	👉👉👉👉👉	2,100 ml	30 pieces	90 pieces (3 x 30 pieces)

medi-inn® CLASSIC Incontinence Slips

Maximum protection – safe all day, every day



- Shaped disposable Incontinence slips
- With wetness indicator
- Unisex fit
- Odour-binding and breathable
- Latex-free
- Fixation: adjustable, resealable tabs on the sides
- Colour coded packaging for easy absorbency level identification
- Non-sterile
- Class I medical device in accordance with Regulation (EU) 2017/745




Item number	Product	Colour coding for absorbency levels	Absorbency in accordance with ISO 11948-1	Multipack	Pieces per carton
N18557	medi-inn® CLASSIC Slips NORMAL M	4 blue droplets	1,950 ml	30 pieces	120 pieces (4 x 30 pieces)
N18558	medi-inn® CLASSIC Slips NORMAL L	4 blue droplets	2,500 ml	30 pieces	120 pieces (4 x 30 pieces)
N18561	medi-inn® CLASSIC Slips PLUS M	5 blue droplets	2,450 ml	30 pieces	90 pieces (3 x 30 pieces)
N18562	medi-inn® CLASSIC Slips PLUS L	5 blue droplets	3,100 ml	15 pieces	60 pieces (4 x 15 pieces)
N18563	medi-inn® CLASSIC Slips PLUS XL	5 blue droplets	3,150 ml	15 pieces	60 pieces (4 x 15 pieces)
N18564	medi-inn® CLASSIC Slips SUPER S	6 green droplets	2,000 ml	30 pieces	120 pieces (4 x 30 pieces)
N18565	medi-inn® CLASSIC Slips SUPER M	6 green droplets	3,150 ml	15 pieces	60 pieces (4 x 15 pieces)
N18566	medi-inn® CLASSIC Slips SUPER L	6 green droplets	3,700 ml	15 pieces	60 pieces (4 x 15 pieces)

medi-inn® CLASSIC Pants



Reliable protection with the comfort of underwear

- Disposable incontinence pants with high absorption capacity
- High elasticated waistband for a comfortable fit
- Unisex fit - can be worn like normal underwear
- Odour-binding and breathable
- Latex-free
- Colour coded packaging for easy absorbency level identification
- FSC®-certified materials
- Dermatologically tested – Dermatest® seal "EXCELLENT"
- Non-sterile
- Class I medical device in accordance with Regulation (EU) 2017/745
- OEKO-TEX® Standard 100



Item number	Product	Colour coding for absorbency levels	Absorbency in accordance with ISO 11948-1	Hip circumference	Multipack	Pieces per carton
N18249	medi-inn® CLASSIC Pants PLUS M	👉👉👉👉👉	1,500 ml	110 - 140 cm	14 pieces	56 pieces (4 x 14 pieces)
N18250	medi-inn® CLASSIC Pants PLUS L	👉👉👉👉👉	1,500 ml	120 - 160 cm	14 pieces	56 pieces (4 x 14 pieces)
N18251	medi-inn® CLASSIC Pants PLUS XL	👉👉👉👉👉	1,500 ml	120 - 160 cm	14 pieces	56 pieces (4 x 14 pieces)
N18245	medi-inn® CLASSIC Pants SUPER S	👈👈👈👈👈👈	1,900 ml	65 - 90 cm	14 pieces	56 pieces (4 x 14 pieces)
N18246	medi-inn® CLASSIC Pants SUPER M	👈👈👈👈👈👈	1,900 ml	80 - 120 cm	14 pieces	56 pieces (4 x 14 pieces)
N18247	medi-inn® CLASSIC Pants SUPER L	👈👈👈👈👈👈	1,900 ml	110 - 140 cm	14 pieces	56 pieces (4 x 14 pieces)
N18248	medi-inn® CLASSIC Pants SUPER XL	👈👈👈👈👈👈	1,900 ml	120 - 160 cm	14 pieces	56 pieces (4 x 14 pieces)

medi-inn® CLASSIC Fix Soft Fixation Pants



Secure support for optimum protection

- Reusable elasticated fixation pants
- For secure fixation of absorbent incontinence products such as medi-inn® CLASSIC shaped pads
- Washable at 60°C in a laundry bag
- Unisex fit
- Latex-free
- Colour coded on the inside for easy size differentiation
- Non-sterile
- Class I medical device in accordance with Regulation (EU) 2017/745
- OEKO-TEX® Standard 100



Item number	Product	Size	Hip circumference	Multipack	Pieces per carton
N18579	medi-inn® CLASSIC Fix Soft	M ●	70 - 110 cm	25 pieces	100 pieces
N18580	medi-inn® CLASSIC Fix Soft	L ●	75 - 120 cm	25 pieces	100 pieces
N18581	medi-inn® CLASSIC Fix Soft	XL ●	80 - 130 cm	25 pieces	100 pieces
N18582	medi-inn® CLASSIC Fix Soft	XXL ●	85 - 140 cm	25 pieces	100 pieces



medi-inn® CLASSIC Underpads SAP Plus Premium

Fully sealed disposable underpads with an absorbent core made of cellulose fluff and superabsorbent polymer (SAP), featuring a liquid-proof and non-slip backing film

- Absorbent core made from cellulose fluff & SAP binds liquids and unpleasant odours
- High absorption capacity
- Liquid-proof, non-slip backing for secure positioning in bed
- Tear-resistant and quick-drying surface
- Materials:
 - Top layer: Polypropylene
 - Back sheet: Polyethylene
 - Absorbent core: Cellulose fluff & SAP
 - Adhesive: hot melt adhesive
- Non-sterile
- Class I medical device in accordance with Regulation (EU) 2017/745



Item number	Absorbency in accordance with ISO 11948-1	Weight	Colour	Dimensions	Multipack	Pieces per carton	Cartons per pallet
N11069	475 ml	30 g	Blue	Length: 60 ± 2 mm Width: 40 ± 2 mm	15 pieces	300 pieces	12
N11071	600 ml	40 g	Blue	Length: 60 ± 2 mm Width: 60 ± 2 mm	15 pieces	150 pieces	18
N11224	860 ml	58 g	Blue	Length: 90 ± 2 mm Width: 60 ± 4 mm	25 pieces	100 pieces	21
730224	1,375 ml	103 g	Blue	Length: 175 ± 2 mm Width: 70 ± 2 mm	15 pieces	90 pieces	16



medi-inn® PRO Underpads PRO DRY

Highly absorbent and breathable disposable underpads with an absorbent core made of cellulose fluff and superabsorbent polymer (SAP), featuring a liquid-proof and non-slip backing film

- Absorbent core made from cellulose fluff & SAP binds liquids and unpleasant odours
- Very high absorption capacity
- Liquid-proof, non-slip backing for secure positioning in bed
- Breathable top layer enhances patient comfort
- Materials:
 - Top layer: Polypropylene
 - Back sheet: Polyethylene
 - Absorbent core: Cellulose fluff & SAP
 - Adhesive: hot melt adhesive
- Non-sterile
- Class I medical device in accordance with Regulation (EU) 2017/745




Item number	Absorbency in accordance with ISO 11948-1	Weight	Colour	Dimensions	Multipack	Pieces per carton	Cartons per pallet
N11983	1,800 ml	100 g	White	Length: 90 ± 2 mm Width: 60 ± 3 mm	10 pieces	100 pieces	12



medi-inn® CLASSIC Bed Protection Pads, With or Without Wings, 5-layer

Washable bed protection pads provide reliable protection for both patient and mattress

- Quick-drying, brushed top layer made of high-performance fibers
- Liquid-proof backing
- Reusable, washable up to 250 times at 95°C
- Suitable for tumble dryers
- Available with or without wings
- Materials: 80% polyester, 20% polyurethane
- Non-sterile
- Class I medical device in accordance with Regulation (EU) 2017/745
- OEKO-TEX® Standard 100




Item number	Absorption capacity	Liquid absorption in accordance with MDS 12/2015	Layers	Design	Colour	Dimensions	Single pack	Pieces per carton	Cartons per pallet
N13198-g	> 3,000 ml / m ²	2,830 g	5-layer	with wings	Green	85 x 90 cm	1 piece	30 pieces	6
N13198-w	> 3,000 ml / m ²	2,830 g	5-layer	with wings	White	85 x 90 cm	1 piece	30 pieces	6
N13199-g	> 3,000 ml / m ²	2,830 g	5-layer	without wings	Green	85 x 90 cm	1 piece	30 pieces	6
N13199-w	> 3,000 ml / m ²	2,830 g	5-layer	without wings	White	85 x 90 cm	1 piece	30 pieces	6




medi-inn® CLASSIC Bed Protection Pads, With or Without Wings, 3 or 4-layer

Washable bed pads offer reliable protection for the patient and the mattress

- Quick-drying, brushed top layer made of high-performance fibers
- Liquid-proof backing
- Reusable, washable up to 250 times at 95°C
- Suitable for tumble dryers
- Available with or without wings
- Materials:
 - 3-layer: 70% polyester, 30% polyurethane
 - 4-layer: 75% polyester, 25% polyurethane
- Non-sterile
- Class I medical device in accordance with Regulation (EU) 2017/745
- OEKO-TEX® Standard 100



Item number	Absorption capacity	Liquid absorption in accordance with MDS 12/2015	Layers	Design	Colour	Dimensions	Single pack	Pieces per carton	Cartons per pallet
N14344-hb	> 2,600 ml / m ²	2,342 g	4-layer	with wings	Light blue	85 x 90 cm	1 piece	30 pieces	6
N14345-hb	> 2,600 ml / m ²	2,342 g	4-layer	without wings	Light blue	85 x 90 cm	1 piece	30 pieces	8
N14994-w	> 2,000 ml / m ²	1,877 g	3-layer	without wings	White	85 x 90 cm	1 piece	30 pieces	12
N14995-w	> 2,000 ml / m ²	1,877 g	3-layer	with wings	White	85 x 90 cm	1 piece	30 pieces	8



Urology

- Urine Bags
- Urine Measurement System
- Silicone Balloon Catheter



medi-inn® PRO Urine bags

Closed urinary drainage system

- Closed urinary drainage system with needle-free urine sampling port
- Universal stepped connector ensures secure catheter connection
- Needle-free sampling port reduces the risk of needlestick injuries and allows easy cleaning due to its smooth design
- Large drip chamber with integrated anti-reflux valve and constant pressure equalisation through the hydrophobic-coated microfilter reduces the risk of urinary tract infections
- Sturdy housing (“doghouse”) ensures secure fixation of the drainage outlet
- Hygienic, convenient one-handed operation thanks to urine drainage with cross valve
- Secure fixation with flexible, stable universal double hanger and additional strap with hanging loop
- Kink-resistant, flexible tubing with clamp and sheet clip
- Sterile
- Free from phthalates
- Cannot be re-sterilised
- Latex-free
- Class Is medical product sterile in accordance with Regulation (EU) 2017/745



CE 2797 MD

Item number	Volume	Tube length	Multipack	Pieces per carton	Cartons per pallet
N19008	2000 ml	120 cm	25 pieces	25 pieces	24
N18983	2000 ml	90 cm	25 pieces	25 pieces	24

medi-inn® PRO Urine Measurement System



Closed single-chamber urine measurement system

- Closed single-chamber urine measuring system with fine graduation on the urimeter for precise monitoring and measurement of the urine volume
- Measuring accuracy urimeter: ± 1 ml with additional drainage outlet on the urimeter
- Universal stepped connector ensures secure catheter connection
- Large drip chamber with integrated anti-reflux valve and constant pressure equalisation through the hydrophobic-coated microfilter reduces the risk of urinary tract infections
- Sturdy housing ("doghouse") ensures secure fixation of the drainage outlet
- Needle-free sampling port reduces the risk of needlestick injuries and allows easy cleaning due to its smooth design
- Large measurement scale with clearly readable markings
- Secure fixation with flexible, stable universal double hanger and additional strap with hanging loop
- Kink-resistant, flexible tubing with locking clamp and sheet clip
- Anti-kink spiral at the proximal tube end reduces the risk of kinking
- Sterile
- Free from phthalates
- Cannot be re-sterilised
- Latex-free
- Class Is medical product sterile in accordance with Regulation (EU) 2017/745




Item number	Volume	Tube length	Multipack	Pieces per carton	Cartons per pallet
N18984	Bag: 2600 ml Measuring chamber: 500 ml	150 cm	10 pieces	10 pieces	20

medi-inn® PRO Silicone Balloon Catheter, 2-Way Tiemann

2-way Tiemann catheter with two functional channels and a curved Tiemann tip – ideal for safe catheterisation in male patients with urinary outflow obstruction

- Made of 100% biocompatible silicone - particularly soft, flexible and well tolerated
- Tiemann tip – specially moulded for easier passage for prostatic hyperplasia or urethral stricture
- Balloon for secure placement in the bladder, easily inflated via plastic valve
- Colour coding for quick and reliable size identification
- Radiopaque contrast line for precise position control in the imaging process
- Individually packed
- Sterile
- Cannot be re-sterilised
- Class IIb medical device sterile in accordance with Regulation (EU) 2017/745




CE 2797 MD   

Item number	Balloon volume in ml	Catheter size (CH / FR)	Outer diameter	Length	Colour code	Multipack	Pieces per carton
N640597	10 ml	12 CH / FR	4.0 mm	350 mm	White	10 pieces	500 pieces
N640598	10 ml	14 CH / FR	4.7 mm	350 mm	Green	10 pieces	500 pieces
N640599	10 ml	16 CH / FR	5.3 mm	350 mm	Orange	10 pieces	500 pieces
N640600	10 ml	18 CH / FR	6.0 mm	350 mm	Red	10 pieces	500 pieces
N640601	10 ml	20 CH / FR	6.7 mm	350 mm	Yellow	10 pieces	500 pieces

medi-inn® PRO Silicone Balloon Catheter, 2-Way Standard

2-way balloon catheter with two functional channels and a straight standard tip – ideal for safe catheterisation in male and female patients with normal urethral anatomy

- Made of 100% biocompatible silicone - particularly soft, flexible and well tolerated
- Rounded standard tip – facilitates atraumatic insertion in normal anatomy
- Balloon for secure placement in the bladder, easily inflated via plastic valve
- Colour coding for quick and reliable size identification
- Radiopaque contrast line for precise position control in the imaging process
- Individually packed
- Sterile
- Cannot be re-sterilised
- Class IIb medical device sterile in accordance with Regulation (EU) 2017/745



Item number	Balloon volume in ml	Catheter size (CH / FR)	Outer diameter	Length	Colour code	Multipack	Pieces per carton
N640602	15 ml	12 CH / FR	4.0 mm	334 mm	White	10 pieces	500 pieces
N640603	15 ml	14 CH / FR	4.7 mm	334 mm	Green	10 pieces	500 pieces
N640604	15 ml	16 CH / FR	5.3 mm	334 mm	Orange	10 pieces	500 pieces
N640605	30 ml	18 CH / FR	6.0 mm	334 mm	Red	10 pieces	500 pieces
N640606	30 ml	20 CH / FR	6.7 mm	334 mm	Yellow	10 pieces	500 pieces
N640607	30 ml	22 CH / FR	7.3 mm	334 mm	Purple	10 pieces	500 pieces
N640608	30 ml	24 CH / FR	8.0 mm	334 mm	Blue	10 pieces	500 pieces

medi-inn® PRO Silicone Balloon Catheter, 3-Way Standard

3-way balloon catheter with three functional channels and straight standard tip – ideal for urinary drainage with additional irrigation channel, e.g. after urological procedures or in cases of haematuria

- The additional irrigation channel enables targeted bladder irrigation after surgical interventions or in cases of haematuria
- Made of 100% biocompatible silicone – particularly soft, flexible and well tolerated
- Balloon for secure placement in the bladder, easily inflated via integrated plastic valve
- Colour-coded size coding for quick identification
- Radiopaque contrast line for precise position control in the imaging process
- Individually packed
- Sterile
- Cannot be re-sterilised
- Class IIb medical device sterile in accordance with Regulation (EU) 2017/745



CE 2797 MD   

Item number	Balloon volume in ml	Catheter size (CH / FR)	Outer diameter	Length	Colour code	Multipack	Pieces per carton
N640611	30 ml	16 CH / FR	5.3 mm	334 mm	Orange	10 pieces	500 pieces
N640612	30 ml	18 CH / FR	6.0 mm	350 mm	Red	10 pieces	500 pieces
N640613	30 ml	20 CH / FR	6.7 mm	350 mm	Yellow	10 pieces	500 pieces
N640614	30 ml	22 CH / FR	7.3 mm	350 mm	Purple	10 pieces	500 pieces
N640615	30 ml	24 CH / FR	8.0 mm	350 mm	Blue	10 pieces	500 pieces

medi-inn® PRO Silicone Balloon Catheter for Children


2-way balloon catheter with two functional channels and a straight standard tip – specially developed for the safe and gentle catheterisation in paediatric patients with normal urethral anatomy

- Made of 100% biocompatible silicone - particularly soft, flexible and well suited for pediatric anatomy
- Smooth surface and atraumatic tip to minimise mechanical irritation
- Balloon for secure placement in the bladder, easily inflated via integrated plastic valve
- Colour-coded size coding for quick identification
- Radiopaque contrast line for precise position control in the imaging process
- Individually packed
- Sterile
- Cannot be re-sterilised
- Class IIb medical device sterile in accordance with Regulation (EU) 2017/745



CE 2797 MD   

Item number	Balloon volume in ml	Catheter size (CH / FR)	Outer diameter	Length	Colour code	Multipack	Pieces per carton
N640616	1.5 ml	6 CH / FR	2.0 mm	234 mm	Pink	10 pieces	500 pieces
N640617	3-5 ml	8 CH / FR	2.7 mm	234 mm	Black	10 pieces	500 pieces
N640618	3-5 ml	10 CH / FR	3.3 mm	234 mm	Grey	10 pieces	500 pieces



10

Surgery

- Sterile Extraction Bags with Handle
- Sterile Extraction Bags without Handle



medi-inn® PRO Specimen Retrieval Bag with Handle

Self-opening, disposable specimen retrieval bag with integrated handle

- Self-opening and releasable specimen retrieval bag with integrated handle
- Transparent, highly tear-resistant, and liquid-tight retrieval bag made of thermoplastic polyurethane (TPU) for the protected collection of small and large resected tissue
- Kink-resistant Nitinol wire for safe intraabdominal closure by pulling the integrated thread
- Colour-coded Medi-Bag handles for quick orientation when choosing the right bag volume from 200 to 1500 ml
- Extremely easy to use thanks to infinitely variable intraoperative opening and closing
- 10 mm trocars for simple insertion of the retrieval bag
- Latex-free
- Single-use
- Double sterile packaging
- Class IIa medical device in accordance with Directive 93/42/EEC


Areas of application

- Adrenalectomy
- Appendectomy
- Cholecystectomy
- Cystectomy
- Ectopic pregnancy
- Laparoscopic colon surgery (bowel, sigmoid and colon resection)
- Lymph node removal
- Nephrectomy
- Prostatectomy
- Thoracoscopy (extraction of lung segments)

CE0123



STERILE



Item number	Materials	Colour	Retrieval bag volume	Retrieval bag diameter	Trocar diameter	SU	Pieces per carton	Cartons per pallet
BP-40350006	Bag: polyurethane (TPU) Closure wire: Nitinol	Pink	200 ml	80 mm	10 mm	5	60 pieces	15
BP- 40350007	Bag: polyurethane (TPU) Closure wire: Nitinol	Blue	500 ml	100 mm	10 mm	5	60 pieces	15
BP- 40350008	Bag: polyurethane (TPU) Closure wire: Nitinol	Green	800 ml	130 mm	10 mm	5	60 pieces	15
BP- 40350010	Bag: polyurethane (TPU) Closure wire: Nitinol	Orange	1500 ml	130 mm	10 mm	5	60 pieces	15

medi-inn® PRO Specimen Retrieval Bag

Self-opening specimen retrieval bag with insertion sheath and pusher. Without handle

- Transparent, highly tear-resistant, and liquid-tight retrieval bag made of thermoplastic polyurethane (TPU) for the protected collection of small and large resected tissue
- Coated and kink-resistant Nitinol wire with memory effect for safe intra-abdominal closure
- 10 and 12 mm trocars for simple insertion of the retrieval bag
- Latex-free
- Single-use
- Double sterile packaging
- Class IIa medical device in accordance with Directive 93/42/EEC

Areas of application

- Adrenalectomy
- Appendectomy
- Cholecystectomy
- Cystectomy
- Ectopic pregnancy
- Laparoscopic colon surgery (bowel, sigmoid and colon resection)
- Lymph node removal
- Nephrectomy
- Prostatectomy
- Thoracoscopy (extraction of lung segments)



CE0123



STERILE EO

Item number	Materials	Colour	Retrieval bag volume	Retrieval bag diameter	Trocar diameter	SU	Pieces per carton	Cartons per pallet
BP-40350011	Bag: polyurethane (TPU) Closure wire: Nitinol	Pink	200 ml	50 mm	10 mm	10	120 pieces	12
BP-40350012	Bag: polyurethane (TPU) Closure wire: Nitinol	Yellow	400 ml	90 mm	10 mm	10	120 pieces	12
BP-40350013	Bag: polyurethane (TPU) Closure wire: Nitinol	Green	800 ml	100 mm	10 mm	10	120 pieces	12



Operating Room Supplies

- Preoperative Preparation and Care
- Patient Positioning and Fixation
- Operating Table Covers
- OR Floor Mats
- Surgical Apparel



ME Medical™ Professional Surgical Clipper with Charging Station

Professional surgical clipper with rechargeable battery

- Professional clipper for the removal of body and head hair
- Set consisting of: SC004 razor and SC010 charger
- Suitable for wet and dry shaving
- Twistable shaving head: Ensures safe handling even in hard-to-reach areas
- With a universal or neuro blade for safe removal of body hair without cuts or micro-injuries
- Non-slip grip and ergonomic design for secure handling
- On/off switch with illuminated display
- Advanced, rechargeable lithium battery technology
- Digital LED battery status display
- Battery overcharge protection
- A high battery performance of 150 minutes and short charging time
- Low noise level
- Easy to clean and disinfect
- Waterproof according to IPX-7 standard
- Class I medical device in accordance with Regulation (EU) 2017/745



Item number	Materials	Colour	Dimensions	Single pack	Pieces per carton	Cartons per pallet
SC014	Clipper with charging station Pollutant-free ABS (acrylonitrile-butadiene styrene copolymer) Blade: Stainless steel	White / Green	Clipper with shaving head 158 mm x 46.5 mm x 40.5 mm Charging station: height 68 mm, width 68 mm, depth 74 mm	1 piece	12 pieces	22



ME Medical™ Professional Surgical Clipper

Professional surgical clipper with rechargeable battery

- Professional clipper for the removal of body and head hair
- Suitable for wet and dry shaving
- Twistable shaving head: Ensures safe handling even in hard-to-reach areas
- With a universal or neuro blade for safe removal of body hair without cuts or micro-injuries
- Non-slip grip and ergonomic design for secure handling
- On/off switch with illuminated display
- Advanced, rechargeable lithium battery technology
- Digital LED battery status display
- Battery overcharge protection
- A high battery performance of 150 minutes and short charging time
- Low noise level
- Easy to clean and disinfect
- Waterproof according to IPX-7 standard
- Class I medical device in accordance with Regulation (EU) 2017/745



Item number	Materials	Colour	Dimensions	Single pack	Pieces per carton	Cartons per pallet
SC004	Clipper: Pollutant-free ABS (acrylonitrile butadiene styrene copolymer) Blade: Stainless steel	White / Green	Clipper with shaving head 158 mm x 46.5 mm x 40.5 mm	1 piece	24 pieces	30



ME Medical™ Charging Station

Charging station for ME Medical™ Surgical Clipper



- Charging station with USB-C port
- Can also be wall-mounted
- Easy to clean and disinfect




Item number	Materials	Colour	Dimensions	Single pack	Pieces per carton	Cartons per pallet
SC010	Charging station: Pollutant-free ABS (acrylonitrile butadiene styrene copolymer)	White / Green	Charging station: height 68 mm, width 68 mm, depth 74 mm	1 piece	24 pieces	30



ME Medical™ Shaving Head for Professional Surgical Clipper

ME Medical™ Universal disposable blade

- The universal blade is designed to meet most hair removal and cutting requirements
- Minimises the risk of micro-injuries, infections, skin irritations
- Efficient preoperative preparation, saving time with easy and safe handling
- Safe removal of body hair without cuts and skin injuries
- Clean and clearly visible incision site, supporting proper adhesion of surgical drapes
- Significant reduction in infection rates
- Individually packaged in peel packaging
- Safe removal without risk of injury
- Non-sterile
- Class I medical device in accordance with Regulation (EU) 2017/745




Item number	Colour	Dimensions	Multipack	Pieces per carton	Cartons per pallet
SC002	Light green	Overall dimensions: 46.5 mm x 30.5 mm / 9 mm Cutting width: 40 mm	50 pieces	300 pieces	36



ME Medical™ Shaving Head for Professional Surgical Clipper

ME Medical™ Neuro disposable blade

- The neuro blade is designed for coarse or thick hair in areas such as the scalp or chest
- Minimises the risk of micro-injuries, infections, skin irritations
- Efficient preoperative preparation, saving time with easy and safe handling
- Safe removal of body hair without cuts and skin injuries
- Clean and clearly visible incision site, supporting proper adhesion of surgical drapes
- Significant reduction in infection rates
- Individually packaged in peel packaging
- Safe removal without risk of injury
- Non-sterile
- Class I medical device in accordance with Regulation (EU) 2017/745




Item number	Colour	Dimensions	Multipack	Pieces per carton	Cartons per pallet
SC003	Dark green	Overall dimensions: 46.5 mm x 30.5 mm / 9 mm Cutting width: 40 mm	25 pieces	150 pieces	72



ME Medical™ Shaving Head for Professional Surgical Clipper

ME Medical™ Precision disposable blade

- The precision blade is designed for sensitive areas of the body such as the genital area
- Minimises the risk of micro-injuries, infections, skin irritations
- Efficient preoperative preparation, saving time with easy and safe handling
- Safe removal of body hair without cuts and skin injuries
- Clean and clearly visible incision site, supporting proper adhesion of surgical drapes
- Significant reduction in infection rates
- Individually packaged in peel packaging
- Safe removal without risk of injury
- Non-sterile
- Class I medical device in accordance with Regulation (EU) 2017/745



Item number	Colour	Dimensions	Multipack	Pieces per carton	Cartons per pallet
SC008	Mint	Overall dimensions: 46.5 mm x 37 mm / 9 mm Cutting width: 25 mm	25 pieces	150 pieces	72



medi-inn® CLASSIC Hand Scrub Brush with Nail Cleaner, Sterile

Sterile scrub brush with nail cleaner

- Hygienic cleaning of palms and fingernails
- Soft sponge and bristles
- Suitable for preoperative preparation and general cleansing
- Sterile
- Not re-sterilizable



Item number	Materials	Colour	Dimensions	Multipack	Pieces per carton	Cartons per pallet
N18854	Brush material: Low-density polyethylene (PE) Sponge material: polyurethane(PU) Nail cleaner material: polypropylene (PP)	White	Brush: length 83 mm width 50 mm, height 39 mm Sponge: length 67 mm width 9.0 mm, height 1.8 mm	312 pieces (12 x 26 pieces)	312 pieces	8



medi-inn® CLASSIC Patient Restraint

Patient restraints for securing extremities to protect patients from injury or harm

- Secure fixation of patients wrists and ankles
- Patient restraint with 250 cm fixation strap
- Blue foam with a white for cushioning
- Velcro fastener
- Stable support thanks to metal ring
- Non-sterile
- Packed in pairs in a polybag
- Class I medical device in accordance with Regulation (EU) 2017/745



Item number	Materials	Colour	Dimensions	Single pack	Pieces per carton	Cartons per pallet
N13147	Foam cuff: polyurethane Fastening ring: steel 304 Retaining straps: nylon	Blue / White	Width: 8 cm Length: 36 cm	1 pair	24 pairs	42



medi-inn® CLASSIC Patient Restraint

Patient restraints for securing extremities to protect patients from injury or harm

- Secure fixation of patients wrists and ankles
- Width: 5 cm
- Available in two lengths: 30 cm or 40 cm
- Inner wool flocking provides cushioning
- Velcro fastener
- Stable support thanks to metal ring
- Non-sterile
- Individually packed in a polybag
- Class I medical device in accordance with Regulation (EU) 2017/745




Item number	Materials	Colour	Dimensions	Single pack	Pieces per carton	Cartons per pallet
N10337-w	Polyester acrylic, wool flocked inside, triangular metal ring	White	Width: 5 cm Length: 30 cm	1 piece	200 pieces	24
N10339-w	Polyester acrylic, wool flocked inside, triangular metal ring	White	Width: 5 cm Length: 40 cm	1 piece	200 pieces	26

medi-inn® CLASSIC Operating Table Cover



Liquid-proof and tear-resistant cover with an absorbent nonwoven top layer



- Applications: operating tables, examination tables, patient beds in hospitals and medical practices
- Reliable protection of operating tables and examination tables against contamination
- High fluid absorption capacity
- Liquid-proof and hygienic
- Reduces the risk of cross-contamination for patients and staff
- Individually vacuum packed
- Non-sterile
- Class I medical device in accordance with Regulation (EU) 2017/745




Item number	Materials	Colour	Dimensions	Absorbency*	Multipack	Pieces per carton	Cartons per pallet
N12583-b	Top layer: polypropylene Back sheet: polyethylene Absorbent core: cellulose fluff & SAP Adhesive: hot melt adhesive	Blue	Length: 152 cm ± 2 Width: 100 cm ± 2	2500 ml	10 pieces	50 pieces	15
N12254-b	Top layer: polypropylene Back sheet: polyethylene Absorbent core: cellulose fluff & SAP Adhesive: hot melt adhesive	Blue	Length: 225 cm ± 2 Width: 100 cm ± 2	4050 ml	10 pieces	20 pieces	18

* in accordance with ISO 11948-1



medi-inn® CLASSIC Operating Table Cover with Carry Sheet

Liquid-proof and tear-resistant cover with an absorbent nonwoven top layer



- Applications: operating tables, patient beds in hospitals and medical practices
- Reliable protection of operating tables and examination tables against contamination
- High fluid absorption capacity
- Liquid-proof and hygienic
- Reduces the risk of cross-contamination for patients and staff
- Individually vacuum packed
- Non-sterile
- Class I medical device in accordance with Regulation (EU) 2017/745



Item number	Materials	Colour	Load capacity Carry sheet	Dimensions	Absorbency*	Multipack	Pieces per carton	Cartons per pallet
N12351-b	Top layer: polypropylene Back sheet: polyethylene Absorbent core: cellulose fluff & SAP Adhesive: hot melt adhesive Carry sheet: 100% polyester	Blue	150 kg	Length: 225 cm ± 2 Width: 100 cm ± 2 Carry sheet: 100 cm x 150 cm	4050 ml	50 pieces	50 pieces	4

* in accordance with ISO 11948-1



medi-inn® CLASSIC Operating Table Cover with Anti-Slip Film

Liquid-proof and tear-resistant cover with an absorbent nonwoven top layer



- Applications: operating tables, examination tables, patient beds in hospitals and medical practices
- Reliable protection of operating tables and examination tables against contamination
- Non-slip, suitable for procedures on inclined operating tables
- High fluid absorption capacity
- Liquid-proof and hygienic
- Reduces the risk of cross-contamination for patients and staff
- Individually vacuum packed
- Non-sterile
- Class I medical device in accordance with Regulation (EU) 2017/745




Item number	Materials	Colour	Dimensions	Absorbency*	Multipack	Pieces per carton	Cartons per pallet
N12457-b	Top layer: polypropylene Back sheet: polyethylene Absorbent core: cellulose fluff & SAP Adhesive: hot melt adhesive	Blue	Length: 225 cm ± 2 Width: 100 cm ± 2	4050 ml	10 pieces	20 pieces	18

* in accordance with ISO 11948-1



medi-inn® CLASSIC Floor Mat

Absorbent floor mat for the operating room area, with non-slip backing, non-sterile



- Reduces the risk of slipping in the operating room area
- High fluid absorption capacity
- Non-slip thanks to anti-slip film on both sides
- With liquid-proof backing
- Helps prevent the spread of fluids within the room
- Contributes to reducing the risk of injuries caused by slipping and falling on wet floors, thereby also lowering the potential risk of infection
- Reduces the risk of cross-contamination for patients and staff
- Non-sterile
- Packed in vacuum bags of 10 pieces




Item number	Materials	Colour	Dimensions	Absorbency*	Multipack	Pieces per carton	Cartons per pallet
N12582	Top layer: non-woven polypropylene fabric Absorbent core: Airlaid + SAP Backing: PE film (anti-slip film)	Blue	Length: 121cm ± 3 Width: 80 cm ± 2	2200 ml	10 pieces	50 pieces	18

* In accordance with ISO 11948-1

medi-inn® CLASSIC SMS Jacket

Long-sleeved SMS jackets with press studs, 3 pockets, knitted cuffs and collar



- Lightweight, breathable warming jacket made of SMS non-woven fabric
- Suitable for use in the operating room and on hospital wards
- Fluid-resistant, protects against moisture
- Knitted collar and cuffs, long sleeves
- 5 press studs for easy donning and doffing
- 3 practical pockets for storing utensils
- Non-sterile
- Class I medical device in accordance with Regulation (EU) 2017/745



Item number	Design	Material	Colour	Dimensions	Multipack	Pieces per carton	Cartons per pallet
N640619-S	Long-sleeved	SMS non-woven (three-layer fiber nonwoven) Material weight 40 g/m ²	Blue	Length 80 cm Width: 120 cm	10 pieces	50 pieces	18
N640619-M	Long-sleeved	SMS non-woven (three-layer fiber nonwoven) Material weight 40 g/m ²	Blue	Length 80 cm Width: 130 cm	10 pieces	50 pieces	18
N640619-L	Long-sleeved	SMS non-woven (three-layer fiber nonwoven) Material weight 40 g/m ²	Blue	Length 80 cm Width: 140 cm	10 pieces	50 pieces	18
N640619-XL	Long-sleeved	SMS non-woven (three-layer fiber nonwoven) Material weight 40 g/m ²	Blue	Length 80 cm Width: 152 cm	10 pieces	50 pieces	18
N640619-XXL	Long-sleeved	SMS non-woven (three-layer fiber nonwoven) Material weight 40 g/m ²	Blue	Length 80 cm Width: 170 cm	10 pieces	50 pieces	18



medi-inn® CLASSIC Universal Disposable Socks

Stretchable disposable socks in universal size

- Stretchable socks in universal size – adapt optimally to any foot
- Versatile use
- Lightweight and breathable
- Latex-free
- Without shaped heel
- Non-sterile



Item number	Materials	Colour	Dimensions	Single pack	Pieces per carton	Cartons per pallet
N15278	80% cotton, 20% nylon	White	Width: 44 cm Length: 9 cm	100 pieces	500 pieces	24



12

Dispenser Systems

- Wall-Mounted Dispenser
- Wall-Mounted Dispenser Accessories
- Stainless Steel Dispensers



medi-inn® CLASSIC Wall-Mounted Dispenser with Long Lever

Aluminium wall-mounted dispenser with cover plate, including mounting material

- Durable aluminium dispenser with stainless steel pump
- For soaps, disinfectants and detergents
- Suitable for 500 ml Euro bottles
- Long stainless steel lever
- Ready for installation, including mounting material
- Delivery includes safety cover



Item number	Suitable for	Version	Colour	Materials	Dimensions	Single pack	Pieces per carton	Cartons per pallet
NI0075-STD	500 ml Euro bottles	Long lever	Silver / white	Housing: aluminium Pump: stainless steel / plastic Lever: Stainless steel	Height: 24.2 cm Width: 8.3 cm Depth (including lever): 23.5 cm	1 piece	10 pieces	17

medi-inn® CLASSIC Wall-Mounted Dispenser with Long Lever, Including Empty Bottle

Aluminium wall-mounted dispenser with refillable 500 ml empty bottle, including mounting material



- Durable aluminium dispenser with stainless steel pump
- For soaps, disinfectants and detergents
- Includes 500 ml Euro bottle
- Long stainless steel lever
- Ready for installation, including mounting material
- Delivery includes safety cover



Item number	Suitable for	Version	Colour	Materials	Dimensions	Single pack	Pieces per carton	Cartons per pallet
NI1468-STD	500 ml Euro bottles	Long lever	Silver / white	Housing: aluminium Pump: stainless steel / plastic Lever: Stainless steel	Height: 24.2 cm Width: 8.3 cm Depth (including lever): 22.5 cm	1 piece	10 pieces	17



medi-inn® CLASSIC Wall-Mounted Dispenser with Short Lever, Including Empty Bottle

Aluminium wall-mounted dispenser with refillable 500 ml empty bottle, including mounting material

- Durable aluminium dispenser with stainless steel pump
- For soaps, disinfectants and detergents
- Includes 500 ml Euro bottle
- Short stainless steel lever
- Ready for installation, including mounting material
- Delivery includes safety cover



Item number	Suitable for	Version	Colour	Materials	Dimensions	Single pack	Pieces per carton	Cartons per pallet
N16412-STD	500 ml Euro bottles	Short lever	Silver / white	Housing: aluminium Pump: stainless steel / plastic Lever: Stainless steel	Height: 24.2 cm Width: 8.5 cm Depth (including lever): 16.0 cm	1 piece	10 pieces	20

medi-inn® CLASSIC Wall-Mounted Dispenser with Short Lever

Aluminium wall-mounted dispenser with cover plate, including mounting material



- Durable aluminium dispenser with stainless steel pump
- For soaps, disinfectants and detergents
- Suitable for 1000 ml Euro bottles
- Short stainless steel lever
- Ready for installation, including mounting material
- Delivery includes safety cover



Item number	Suitable for	Version	Colour	Materials	Dimensions	Single pack	Pieces per carton	Cartons per pallet
N12311-STD	1000 ml Euro bottles	Short lever	Silver / white	Housing: aluminium Pump: stainless steel / plastic Lever: Stainless steel	Height: 29.0 cm Width: 9.5 cm Depth (including lever): 16.5 cm	1 piece	10 pieces	16



medi-inn® CLASSIC Wall-Mounted Dispenser with Long Lever

Aluminium wall-mounted dispenser with cover plate, including mounting material

- Durable aluminium dispenser with stainless steel pump
- For soaps, disinfectants and detergents
- Suitable for 1000 ml Euro bottles
- Long stainless steel lever
- Ready for installation, including mounting material
- Delivery includes safety cover



Item number	Suitable for	Version	Colour	Materials	Dimensions	Single pack	Pieces per carton	Cartons per pallet
N16411-STD	1000 ml Euro bottles	Long lever	Silver / white	Housing: aluminium Pump: stainless steel / plastic Lever: Stainless steel	Height: 29.0 cm Width: 9.5 cm Depth (including lever): 23.5 cm	1 piece	10 pieces	12

medi-inn® CLASSIC Wall-Mounted Dispenser with long lever

Closed wall-mounted dispenser made of plastic with viewing window, including mounting material



- Compact plastic dispenser with viewing window and mounting material
- For soaps, disinfectants and detergents
- Suitable for 500 ml Euro bottles
- Long arm lever made of high-quality plastic
- Ready for installation, including mounting material
- Delivery includes safety cover



Item number	Suitable for	Version	Colour	Materials	Dimensions	Single pack	Pieces per carton	Cartons per pallet
780080-STD	500 ml Euro bottles	Long lever	White / grey	Housing: plastic Pump: stainless steel / plastic Lever: plastic	Height: 23.5 cm Width: 8.5 cm Depth (including lever): 21.0 cm	1 piece	10 pieces	21



medi-inn® CLASSIC Wall-Mounted Dispenser with Long Lever

Closed wall-mounted dispenser made of plastic with viewing window, including mounting material

- Compact plastic dispenser with viewing window
- For soaps, disinfectants and detergents
- Suitable for 1000 ml Euro bottles
- Dispenser and long lever made from high-quality plastic
- Ready for installation, including mounting material



Item number	Suitable for	Version	Colour	Materials	Dimensions	Single pack	Pieces per carton	Cartons per pallet
N10826-STD	1000 ml Euro bottles	Long lever	White / grey	Housing: plastic Pump: stainless steel / plastic Lever: plastic	Height: 27.7 cm Width: 10.5 cm	1 piece	10 pieces	16

medi-inn® CLASSIC Wall-Mounted Dispenser with Drip Tray and Long Lever

Closed wall-mounted dispenser made of plastic with viewing window and drip tray, including mounting material



- Compact plastic dispenser with viewing window
- For soaps, disinfectants and washing lotions, with practical drip tray, including mounting material
- Suitable for 500 ml Euro bottles
- Dispenser and long lever made of high-quality plastic
- Ready for installation, including mounting material



Item number	Suitable for	Version	Colour	Materials	Dimensions	Single pack	Pieces per carton	Cartons per pallet
N12310-STD	500 ml Euro bottles	Long lever	White / grey	Housing: plastic Pump: stainless steel / plastic Lever: plastic	Height: 23.0 cm Width: 14.8 cm Depth (including lever): 34.4 cm	1 piece	10 pieces	21



medi-inn® CLASSIC Drip Tray for 500-ml or 1000-ml Wall-Mounted Dispensers

Drip tray for mounting on 500 ml or 1000 ml aluminium wall dispensers

- Stainless steel frame for plastic tray
- For easy mounting above the installed dispenser
- Removable tray for easy cleaning and maintenance
- 10 or 20 pieces per carton



Item number	Suitable for	Colour	Materials	Dimensions	Single pack	Pieces per carton	Cartons per pallet
N13268-STD	500 ml wall dispenser	White / silver	Attachment: stainless steel Tray: plastic	Tray: 18.5 cm x 15.5 cm x 3.6 cm (L x W x H) Weight: 175 g	1 piece	10 pieces	48
N15758	1000 ml wall dispenser	White / silver	Attachment: stainless steel Tray: plastic	Tray: 18.5 cm x 15.5 cm x 3.6 cm (L x W x H) Weight: 175 g	1 piece	20 pieces	32

medi-inn® CLASSIC Replacement Pumps for Aluminium Wall-Mounted Dispensers



Replacement pumps for 500 ml or 1000 ml aluminium wall dispensers

- Stainless steel with plastic tube
- Easy to replace



Item number	Suitable for	Materials	Single pack	Pieces per carton
N13196	500 ml wall dispenser	Stainless steel / plastic	1 piece	100 pieces
N14565	1000 ml wall dispenser	Stainless steel / plastic	1 piece	100 pieces



medi-inn® CLASSIC Replacement Pumps for Plastic Wall-Mounted Dispensers

Replacement pumps for 500 ml or 1000 ml aluminium wall dispensers

- Easy to replace



Item number	Suitable for	Materials	Multipack	Pieces per carton
N10950	500 ml Euro bottles	Stainless steel / plastic	50 pieces	250 pieces
N12330	1000 ml Euro bottles	Stainless steel / plastic	50 pieces	250 pieces

medi-inn® CLASSIC Stainless Steel Dispensers for Disposable Gloves



Stainless steel dispensers for the convenient storage of glove boxes and quick access to disposable gloves

- Dispenser for holding one glove box
- For wall mounting
- Made of stainless steel
- Including mounting material
- Colour: steel grey
- Suitable for standard glove boxes
Size XL suitable for glove box heights between 120 and 135 mm



Item number	Product name	Dimensions Width / height / depth	Capacity	SU	Pieces per carton	Cartons per pallet
N12797-XL	Stainless steel dispenser XL	252 x 135 x 80 mm	1 glove box	1 piece	9 pieces	24

medi-inn® CLASSIC Multi-Box Stainless Steel Dispenser for Disposable Gloves



Stainless steel dispenser for the convenient storage of glove boxes and quick access to disposable gloves

- Multi-Box dispenser for storing up to three or four glove boxes
- For wall mounting
- Made of stainless steel
- Including mounting material
- Colour: steel grey
- Suitable for standard glove boxes
item no. N13218 suitable for glove box widths between 240 and 255 mm
item no. N14097 suitable for glove box widths between 190 and 220 mm
item no. N13219 suitable for glove box widths between 240 and 255 mm



Item example N13219

Item number	Product name	Dimensions Width / height / depth	Capacity	SU	Pieces per carton	Cartons per pallet
N13218	3-way stainless steel dispenser	255 x 390 x 80 mm	3 glove boxes	1 piece	9 pieces	8
N14097	3-way stainless steel dispenser, narrow	220 x 350 x 70 mm	3 glove boxes	1 piece	9 pieces	8
N13219	4-way stainless steel dispenser	255 x 515 x 80 mm	4 glove boxes	1 piece	9 pieces	6



medi-inn® CLASSIC Stainless Steel Holder for Aprons in Box

Wall-mounted holder for aprons in boxes

- Made of stainless steel
- For wall mounting, including mounting material
- Designed for storage and dispensing of aprons in a box
- Compatible with BODY PRODUCTS medi-inn® CLASSIC PE Apron, PLUS, item N12989 & N13461



Item number	Material	Colour	Dimensions	Single pack	Pieces per carton	Cartons per pallet
N13504	Stainless steel	Silver matt	Height: 110 mm Width: 275 mm Depth: 125 mm	1 piece	6 pieces	27

medi-inn® CLASSIC Stainless Steel Holder for Blocked Aprons

Adjustable wall bracket for blocked aprons



- Made of stainless steel
- For wall mounting
- Designed for storage and dispensing of blocked aprons
- Hook spacing individually adjustable



Item number	Material	Colour	Dimensions	Single pack	Pieces per carton	Cartons per pallet
N13505	Stainless steel	Silver matt	Height: 65 mm Width: 360 mm Depth: 10 mm	1 piece	6 pieces	42



medi-inn® CLASSIC Stainless Steel Holder SINGLE for Face Masks

Stainless steel dispenser for one box containing 50 disposable face masks

- For wall mounting, including mounting material
- Made of stainless steel



Item number	Material	Colour	Dimensions	Single pack	Pieces per carton	Cartons per pallet
N18313	Stainless steel	steel grey	Height: 190 mm Width: 108 mm Depth: 85 mm	1 piece	10 pieces	28



Your contact:

BODY PRODUCTS
Relax Pharma und Kosmetik GmbH

Alfred-Nobel-Strasse 1-3
50226 Frechen

Tel.: +49 2234 91165-0
Fax: +49 2234 91165-23
Email: info@body-products.de

WWW.BODY-PRODUCTS.DE